

# 2019 order form

Use dimension chart to select appropriate size.




Mini and small sizes available Spring 2019.

		CHOOSE	NEW K610 mini	NEW K620 small	K630 medium	K640 large	K650 XL		
* Color	CHOOSE	Red	<input type="checkbox"/>	Red	<input type="checkbox"/>	Red	<input type="checkbox"/>	Red	<input type="checkbox"/>
		Blue	<input type="checkbox"/>	Blue	<input type="checkbox"/>	Blue	<input type="checkbox"/>	Blue	<input type="checkbox"/>
		Gray	<input type="checkbox"/>	Gray	<input type="checkbox"/>	Gray	<input type="checkbox"/>	Gray	<input type="checkbox"/>
		Lime	<input type="checkbox"/>	Lime	<input type="checkbox"/>	Lime	<input type="checkbox"/>	Lime	<input type="checkbox"/>
		Pink	<input type="checkbox"/>	Pink	<input type="checkbox"/>	Pink	<input type="checkbox"/>	Pink	<input type="checkbox"/>
* Upper frame	CHOOSE	Dynamic	N/A	K621 \$1,000 <input type="checkbox"/>	K631 \$1,375 <input type="checkbox"/>	K641 \$1,555 <input type="checkbox"/>	K651 \$2,005 <input type="checkbox"/>		
	CHOOSE	Standard	N/A	K622 \$590 <input type="checkbox"/>	K632 \$815 <input type="checkbox"/>	K642 \$980 <input type="checkbox"/>	K652 \$1,435 <input type="checkbox"/>		
* Base	Standard base	CHOOSE	With odometer	K613 \$860 <input type="checkbox"/>	K623 \$530 <input type="checkbox"/>	K633 \$665 <input type="checkbox"/>	K643 \$830 <input type="checkbox"/>	K653 \$1,230 <input type="checkbox"/>	
		CHOOSE	Without odometer	K614 \$710 <input type="checkbox"/>	K624 \$370 <input type="checkbox"/>	K634 \$500 <input type="checkbox"/>	K644 \$670 <input type="checkbox"/>	K654 \$1,060 <input type="checkbox"/>	
	Utility base (larger wheels)	CHOOSE	With odometer	N/A	K635 \$915 <input type="checkbox"/>	K635 \$915 <input type="checkbox"/>	K645 \$1,080 <input type="checkbox"/>	N/A	
		CHOOSE	Without odometer	N/A	K636 \$750 <input type="checkbox"/>	K636 \$750 <input type="checkbox"/>	K646 \$920 <input type="checkbox"/>	N/A	
	Treadmill/stability base (comes in gray)	CHOOSE	Standard	N/A	K676 \$945 <input type="checkbox"/>	K676 \$945 <input type="checkbox"/>	K676 \$945 <input type="checkbox"/>	N/A	
		CHOOSE	Wide	N/A	K678 \$1,060 <input type="checkbox"/>	K678 \$1,060 <input type="checkbox"/>	K678 \$1,060 <input type="checkbox"/>	N/A	
* Arm supports (pair)	CHOOSE	Small arm prompts	K512 \$475 <input type="checkbox"/>	K512 \$475 <input type="checkbox"/>	K512 \$475 <input type="checkbox"/>	K512 \$475 <input type="checkbox"/>	K512 \$475 <input type="checkbox"/>		
		Large arm prompts	N/A	N/A	K532 \$510 <input type="checkbox"/>	K532 \$510 <input type="checkbox"/>	K532 \$510 <input type="checkbox"/>		
		Arm platforms	K697 \$340 <input type="checkbox"/>	K697 \$340 <input type="checkbox"/>	K697 \$340 <input type="checkbox"/>	K697 \$340 <input type="checkbox"/>	K697 \$340 <input type="checkbox"/>		
		Arm platforms w/ handgrips	K696 \$390 <input type="checkbox"/>	K696 \$390 <input type="checkbox"/>	K696 \$390 <input type="checkbox"/>	K696 \$390 <input type="checkbox"/>	K696 \$390 <input type="checkbox"/>		
Handbrake for utility base (requires arm prompts)			N/A	N/A	K687 \$270 <input type="checkbox"/>	K688 \$295 <input type="checkbox"/>	N/A		
Hand loops (pair)			K579 \$185 <input type="checkbox"/>	K579 \$185 <input type="checkbox"/>	K579 \$185 <input type="checkbox"/>	K579 \$185 <input type="checkbox"/>	K579 \$185 <input type="checkbox"/>		
* Chest prompt	CHOOSE	Small	K514 \$510 <input type="checkbox"/>	K514 \$510 <input type="checkbox"/>	K514 \$510 <input type="checkbox"/>	N/A	N/A		
		Medium	K524 \$550 <input type="checkbox"/>	K524 \$550 <input type="checkbox"/>	K524 \$550 <input type="checkbox"/>	K524 \$550 <input type="checkbox"/>	K524 \$550 <input type="checkbox"/>		
		Large	K534 \$595 <input type="checkbox"/>	K534 \$595 <input type="checkbox"/>	K534 \$595 <input type="checkbox"/>	K534 \$595 <input type="checkbox"/>	K534 \$595 <input type="checkbox"/>		
* Multi-position saddle	CHOOSE	NEW! Small	N/A	K673 \$850 <input type="checkbox"/>	N/A	N/A	N/A		
	CHOOSE	Large	N/A	N/A	K674 \$910 <input type="checkbox"/>	K674 \$910 <input type="checkbox"/>	K674 \$910 <input type="checkbox"/>		
* Pelvic support (includes handholds)	CHOOSE	Small	K556 \$385 <input type="checkbox"/>	K556 \$385 <input type="checkbox"/>	K556 \$385 <input type="checkbox"/>	N/A	N/A		
		Medium	N/A	K557 \$405 <input type="checkbox"/>	K557 \$405 <input type="checkbox"/>	K557 \$405 <input type="checkbox"/>	N/A		
		Large	N/A	N/A	K558 \$420 <input type="checkbox"/>	K558 \$420 <input type="checkbox"/>	K558 \$420 <input type="checkbox"/>		
* Hip positioner (includes handholds)	CHOOSE	Small	K513 \$385 <input type="checkbox"/>	K513 \$385 <input type="checkbox"/>	K513 \$385 <input type="checkbox"/>	K513 \$385 <input type="checkbox"/>	K513 \$385 <input type="checkbox"/>		
		Large	K533 \$420 <input type="checkbox"/>	K533 \$420 <input type="checkbox"/>	K533 \$420 <input type="checkbox"/>	K533 \$420 <input type="checkbox"/>	K533 \$420 <input type="checkbox"/>		
* Hip positioner pad	CHOOSE	Small	K541 \$98 <input type="checkbox"/>	K541 \$98 <input type="checkbox"/>	K541 \$98 <input type="checkbox"/>	K541 \$98 <input type="checkbox"/>	K541 \$98 <input type="checkbox"/>		
		Large	N/A	N/A	K542 \$130 <input type="checkbox"/>	K542 \$130 <input type="checkbox"/>	K542 \$130 <input type="checkbox"/>		
* Thigh prompts (pair)	CHOOSE	Small	K515 \$380 <input type="checkbox"/>	K515 \$380 <input type="checkbox"/>	K515 \$380 <input type="checkbox"/>	K515 \$380 <input type="checkbox"/>	N/A		
		Large	N/A	N/A	K535 \$450 <input type="checkbox"/>	K535 \$450 <input type="checkbox"/>	K535 \$450 <input type="checkbox"/>		
* Ankle prompts for standard base (pair)			K681 \$150 <input type="checkbox"/>	K682 \$160 <input type="checkbox"/>	K683 \$175 <input type="checkbox"/>	K684 \$190 <input type="checkbox"/>	K684 \$190 <input type="checkbox"/>		
* Ankle prompts for utility base (pair)			N/A	N/A	K693 \$175 <input type="checkbox"/>	K694 \$195 <input type="checkbox"/>	N/A		
* Ankle prompts for treadmill/stability base (pair)			N/A	N/A	K664 \$195 <input type="checkbox"/>	K664 \$195 <input type="checkbox"/>	N/A		
* Additional handholds (pair)			K519 \$105 <input type="checkbox"/>	K519 \$105 <input type="checkbox"/>	K519 \$105 <input type="checkbox"/>	K519 \$105 <input type="checkbox"/>	K519 \$105 <input type="checkbox"/>		
* Communication tray			K518 \$360 <input type="checkbox"/>	K698 \$350 <input type="checkbox"/>	K698 \$350 <input type="checkbox"/>	K698 \$350 <input type="checkbox"/>	K698 \$350 <input type="checkbox"/>		
* Attendant guide bar			K589 \$130 <input type="checkbox"/>	K679 \$130 <input type="checkbox"/>	K689 \$130 <input type="checkbox"/>	K689 \$130 <input type="checkbox"/>	N/A		
* Rifton accessories tote			K522 \$100 <input type="checkbox"/>	K522 \$100 <input type="checkbox"/>	K522 \$100 <input type="checkbox"/>	K522 \$100 <input type="checkbox"/>	K522 \$100 <input type="checkbox"/>		
<b>Total</b>			<b>K610 \$ _____</b>	<b>K620 \$ _____</b>	<b>K630 \$ _____</b>	<b>K640 \$ _____</b>	<b>K650 \$ _____</b>		


\* Required

# Pacer dimensions

Key user dimensions (inches)		K610 mini	K620 small	K630 medium	K640 large	K650 XL
Elbow height		15½–20½	18½–27½	24–35	32–47	34–49
 <p><b>Key user dimension: elbow height</b> Measure the vertical distance from the bent <b>elbow to the floor</b> while the user is standing upright. Choose the size that allows for growth.</p>						
Item dimensions (inches)		K610 mini	K620 small	K630 medium	K640 large	K650 XL
Floor to top of arm prompt		15½–20½	19–29	24–35	32–47	34–49
Floor to top of MPS (not tilted)*		N/A	10–20	9½–23½	17–35	19½–37½
Overall width	Standard base:	20½	23	26	28	31½
	Utility base:	N/A	N/A	30	32	N/A
	Narrow treadmill/stability base:	N/A	N/A	35	35	N/A
	Wide treadmill/stability base:	N/A	N/A	40	40	N/A
Overall length	Standard base:	22½	26	30	36	39¾
	Utility base:	N/A	N/A	36	40½	N/A
	Treadmill/stability bases:	N/A	N/A	41	41	N/A
Frame height	Standard & utility bases:	12½	16–21	21½–27	28½–38½	31–41
	Treadmill/stability bases:	N/A	23–35	23½–36½	30–47	N/A
Base height (without upper frame)	Standard & utility bases:	N/A	11	14½	16	18½
	Treadmill/stability bases:	N/A	N/A	21	21	N/A
Frame weight (lbs)	Standard base:	9½	10½	15	15½	22½
	Utility base:	N/A	N/A	18½	19¼	N/A
	Narrow treadmill/stability base:	N/A	N/A	27½	27½	N/A
	Wide treadmill/stability base:	N/A	N/A	28	28	N/A
	Standard upper:	N/A	6	7½	9½	11
	Dynamic upper:	N/A	7	11½	13½	15
Dynamic upper movement	Vertical:	N/A	1¼	2	2	2
	Horizontal:	N/A	1½	2	2	2
Maximum treadmill width	Narrow treadmill/stability base:	N/A	N/A	29	29	N/A
	Wide treadmill/stability base:	N/A	N/A	34	34	N/A
Maximum treadmill height	Treadmill/stability bases:	N/A	N/A	11	11	N/A
Maximum working load (lbs)		50	75	150	200	250

\*This measurement is taken at the lowest frame and MPS height, and the highest frame and MPS height.

## Pelvic support

User dimensions (inches) Indicated by the red arrow below			
<b>Key user dimension:</b> half of thigh circumference measured at the groin.			
K556 small	K557 med.	K558 large	
4–8	6–10	8–14	
			
Product dimensions (inches)			
K556 small	K557 med.	K558 large	
<b>a</b>	6½	8½	12
<b>b</b>	12	16	20
<b>c</b>	7	9	12

## Hip positioner

Dimensions (inches) – without pad		
	K513 small	K533 large
<b>a</b> Crotch length	2½	3½
<b>b</b> Width	1¼	1¾
<b>c</b> Length	9	11