

2022 order form

Use dimension chart to select appropriate size.

| | | | K610 mini | K620 small | K630 medium | K640 large | K650 XL | | |
|--|--|----------------------------|-------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|--------------------------|
| * Color | CHOOSE | Red | <input type="checkbox"/> | Red | <input type="checkbox"/> | Red | <input type="checkbox"/> | Red | <input type="checkbox"/> |
| | | Blue | <input type="checkbox"/> | Blue | <input type="checkbox"/> | Blue | <input type="checkbox"/> | Blue | <input type="checkbox"/> |
| | | Gray | <input type="checkbox"/> | Gray | <input type="checkbox"/> | Gray | <input type="checkbox"/> | Gray | <input type="checkbox"/> |
| | | Lime | <input type="checkbox"/> | Lime | <input type="checkbox"/> | Lime | <input type="checkbox"/> | Lime | <input type="checkbox"/> |
| | | Pink | <input type="checkbox"/> | Pink | <input type="checkbox"/> | Pink | <input type="checkbox"/> | Pink | <input type="checkbox"/> |
| * Upper frame | CHOOSE | Dynamic | N/A | K621 \$1,100 <input type="checkbox"/> | K631 \$1,525 <input type="checkbox"/> | K641 \$1,715 <input type="checkbox"/> | K651 \$2,205 <input type="checkbox"/> | | |
| | | Standard | N/A | K622 \$650 <input type="checkbox"/> | K632 \$895 <input type="checkbox"/> | K642 \$1,080 <input type="checkbox"/> | K652 \$1,575 <input type="checkbox"/> | | |
| * Base | Standard base | CHOOSE | With odometer | K613 \$935 <input type="checkbox"/> | K623 \$585 <input type="checkbox"/> | K633 \$730 <input type="checkbox"/> | K643 \$915 <input type="checkbox"/> | K653 \$1,340 <input type="checkbox"/> | |
| | | | Without odometer | K614 \$780 <input type="checkbox"/> | K624 \$410 <input type="checkbox"/> | K634 \$550 <input type="checkbox"/> | K644 \$740 <input type="checkbox"/> | K654 \$1,165 <input type="checkbox"/> | |
| | Utility base (larger wheels) | CHOOSE | With odometer | N/A | N/A | K635 \$1,010 <input type="checkbox"/> | K645 \$1,185 <input type="checkbox"/> | N/A | |
| | | | Without odometer | N/A | N/A | K636 \$830 <input type="checkbox"/> | K646 \$1,015 <input type="checkbox"/> | N/A | |
| | Treadmill/stability base (comes in gray only) | CHOOSE | Standard | N/A | K676 \$1,040 <input type="checkbox"/> | K676 \$1,040 <input type="checkbox"/> | K676 \$1,040 <input type="checkbox"/> | N/A | |
| | | | Wide | N/A | K678 \$1,165 <input type="checkbox"/> | K678 \$1,165 <input type="checkbox"/> | K678 \$1,165 <input type="checkbox"/> | N/A | |
| Arm supports (pair) | CHOOSE | Small arm prompts | K512 \$510 <input type="checkbox"/> | K512 \$510 <input type="checkbox"/> | K512 \$510 <input type="checkbox"/> | K512 \$510 <input type="checkbox"/> | K512 \$510 <input type="checkbox"/> | | |
| | | Large arm prompts | N/A | N/A | K532 \$545 <input type="checkbox"/> | K532 \$545 <input type="checkbox"/> | K532 \$545 <input type="checkbox"/> | | |
| | | Arm platforms | N/A | N/A | K697 \$370 <input type="checkbox"/> | K697 \$370 <input type="checkbox"/> | K697 \$370 <input type="checkbox"/> | | |
| | | Arm platforms w/ handgrips | N/A | N/A | K696 \$420 <input type="checkbox"/> | K696 \$420 <input type="checkbox"/> | K696 \$420 <input type="checkbox"/> | | |
| Handbrake for utility base (requires arm supports) | | | N/A | N/A | K687 \$300 <input type="checkbox"/> | K688 \$325 <input type="checkbox"/> | N/A | | |
| Handlebars (pair) | CHOOSE | Small | K625 \$195 <input type="checkbox"/> | K625 \$195 <input type="checkbox"/> | K625 \$195 <input type="checkbox"/> | K625 \$195 <input type="checkbox"/> | K625 \$195 <input type="checkbox"/> | | |
| | | Large | N/A | N/A | K675 \$205 <input type="checkbox"/> | K675 \$205 <input type="checkbox"/> | K675 \$205 <input type="checkbox"/> | | |
| Hand loops (pair) | | | K579 \$168 <input type="checkbox"/> | K579 \$168 <input type="checkbox"/> | K579 \$168 <input type="checkbox"/> | K579 \$168 <input type="checkbox"/> | K579 \$168 <input type="checkbox"/> | | |
| Chest prompt | CHOOSE | Small | K514 \$535 <input type="checkbox"/> | K514 \$535 <input type="checkbox"/> | K514 \$535 <input type="checkbox"/> | N/A | N/A | | |
| | | Medium | N/A | K524 \$580 <input type="checkbox"/> | K524 \$580 <input type="checkbox"/> | K524 \$580 <input type="checkbox"/> | K524 \$580 <input type="checkbox"/> | | |
| | | Large | N/A | N/A | K534 \$625 <input type="checkbox"/> | K534 \$625 <input type="checkbox"/> | K534 \$625 <input type="checkbox"/> | | |
| Chest pad | | | N/A | K629 \$295 <input type="checkbox"/> | K639 \$315 <input type="checkbox"/> | K649 \$335 <input type="checkbox"/> | N/A | | |
| Multi-Position Saddle | CHOOSE | Small | N/A | K672 \$895 <input type="checkbox"/> | N/A | N/A | N/A | | |
| | | Medium | N/A | N/A | K673 \$920 <input type="checkbox"/> | N/A | N/A | | |
| | | Large | N/A | N/A | K674 \$955 <input type="checkbox"/> | K674 \$955 <input type="checkbox"/> | K674 \$955 <input type="checkbox"/> | | |
| MPS seat cover | CHOOSE | Small/Medium | N/A | K627 \$75 <input type="checkbox"/> | K627 \$75 <input type="checkbox"/> | N/A | N/A | | |
| | | Large | N/A | N/A | N/A | K628 \$75 <input type="checkbox"/> | K628 \$75 <input type="checkbox"/> | | |
| Pelvic support (includes handholds) | CHOOSE | Small | K556 \$415 <input type="checkbox"/> | K556 \$415 <input type="checkbox"/> | K556 \$415 <input type="checkbox"/> | N/A | N/A | | |
| | | Medium | N/A | K557 \$435 <input type="checkbox"/> | K557 \$435 <input type="checkbox"/> | K557 \$435 <input type="checkbox"/> | N/A | | |
| | | Large | N/A | N/A | K558 \$450 <input type="checkbox"/> | K558 \$450 <input type="checkbox"/> | K558 \$450 <input type="checkbox"/> | | |
| Hip positioner (includes handholds) | CHOOSE | Small | K513 \$415 <input type="checkbox"/> | K513 \$415 <input type="checkbox"/> | K513 \$415 <input type="checkbox"/> | K513 \$415 <input type="checkbox"/> | K513 \$415 <input type="checkbox"/> | | |
| | | Large | N/A | N/A | K533 \$450 <input type="checkbox"/> | K533 \$450 <input type="checkbox"/> | K533 \$450 <input type="checkbox"/> | | |
| Hip positioner pad | CHOOSE | Small | K541 \$105 <input type="checkbox"/> | K541 \$105 <input type="checkbox"/> | K541 \$105 <input type="checkbox"/> | K541 \$105 <input type="checkbox"/> | K541 \$105 <input type="checkbox"/> | | |
| | | Large | N/A | N/A | K542 \$142 <input type="checkbox"/> | K542 \$142 <input type="checkbox"/> | K542 \$142 <input type="checkbox"/> | | |
| Thigh prompts (pair) | CHOOSE | Small | N/A | K515 \$400 <input type="checkbox"/> | K515 \$400 <input type="checkbox"/> | K515 \$400 <input type="checkbox"/> | N/A | | |
| | | Large | N/A | N/A | K535 \$470 <input type="checkbox"/> | K535 \$470 <input type="checkbox"/> | K535 \$470 <input type="checkbox"/> | | |
| Ankle prompts for standard base (pair) | | | K681 \$168 <input type="checkbox"/> | K682 \$178 <input type="checkbox"/> | K683 \$195 <input type="checkbox"/> | K684 \$210 <input type="checkbox"/> | K684 \$210 <input type="checkbox"/> | | |
| Ankle prompts for utility base (pair) | | | N/A | N/A | K693 \$195 <input type="checkbox"/> | K694 \$215 <input type="checkbox"/> | N/A | | |
| Ankle prompts for treadmill/stability base (pair) | | | N/A | K664 \$215 <input type="checkbox"/> | K664 \$215 <input type="checkbox"/> | K664 \$215 <input type="checkbox"/> | N/A | | |
| Additional handholds (pair) | | | K519 \$115 <input type="checkbox"/> | K519 \$115 <input type="checkbox"/> | K519 \$115 <input type="checkbox"/> | K519 \$115 <input type="checkbox"/> | K519 \$115 <input type="checkbox"/> | | |
| Communication tray | | | K518 \$385 <input type="checkbox"/> | K698 \$390 <input type="checkbox"/> | K698 \$390 <input type="checkbox"/> | K698 \$390 <input type="checkbox"/> | K698 \$390 <input type="checkbox"/> | | |
| Attendant guide bar | | | K589 \$147 <input type="checkbox"/> | K679 \$147 <input type="checkbox"/> | K689 \$147 <input type="checkbox"/> | K689 \$147 <input type="checkbox"/> | N/A | | |
| Rifton accessories tote | | | K522 \$130 <input type="checkbox"/> | K522 \$130 <input type="checkbox"/> | K522 \$130 <input type="checkbox"/> | K522 \$130 <input type="checkbox"/> | K522 \$130 <input type="checkbox"/> | | |
| Total | | | K610 \$ _____ | K620 \$ _____ | K630 \$ _____ | K640 \$ _____ | K650 \$ _____ | | |


* Required

Pacer dimensions


| Key user dimensions (inches) | | K610 mini | K620 small | K630 medium | K640 large | K650 XL |
|---|---|-----------|------------|-------------|------------|---------|
| Elbow height | | 15½–20½ | 18½–27½ | 24–35 | 32–47 | 34–49 |
|  | Key user dimension: elbow height Measure the vertical distance from the bent elbow to the floor while the user is standing upright. Choose the size that allows for growth. | | | | | |
| Item dimensions (inches) | | K610 mini | K620 small | K630 medium | K640 large | K650 XL |
| Floor to top of arm prompt | | 15½–20½ | 18½–27½ | 24–35 | 32–47 | 34–49 |
| Floor to top of MPS (not tilted)* | | N/A | 10–20 | 9½–23½ | 17–35 | 19½–37½ |
| Overall width | Standard base: | 20½ | 23 | 26 | 28 | 31½ |
| | Utility base: | N/A | N/A | 30 | 32 | N/A |
| | Narrow treadmill/stability base: | N/A | 35 | 35 | 35 | N/A |
| | Wide treadmill/stability base: | N/A | 40 | 40 | 40 | N/A |
| Overall length | Standard base: | 22½ | 26 | 30 | 36 | 39¾ |
| | Utility base: | N/A | N/A | 36 | 40½ | N/A |
| | Treadmill/stability bases: | N/A | 41 | 41 | 41 | N/A |
| Frame height | Standard & utility bases: | 12½ | 16–21 | 21½–27 | 28½–38½ | 31–41 |
| | Treadmill/stability bases: | N/A | 22–34 | 23½–36 | 29½–46½ | N/A |
| Base height (without upper frame) | Standard & utility bases: | N/A | 11 | 14½ | 16 | 18½ |
| | Treadmill/stability bases: | N/A | 21½–28 | 21½–28 | 21½–28 | N/A |
| Frame weight (lb) | Standard base: | 10¼ | 11½ | 15 | 15½ | 22½ |
| | Utility base: | N/A | N/A | 18½ | 19¼ | N/A |
| | Narrow treadmill/stability base: | N/A | 27½ | 27½ | 27½ | N/A |
| | Wide treadmill/stability base: | N/A | 28 | 28 | 28 | N/A |
| | Standard upper: | N/A | 6 | 7½ | 9½ | 11 |
| Dynamic upper movement | Vertical: | N/A | 1¼ | 2 | 2 | 2 |
| | Horizontal: | N/A | 1½ | 2 | 2 | 2 |
| Maximum treadmill width | Narrow treadmill/stability base: | N/A | 29 | 29 | 29 | N/A |
| | Wide treadmill/stability base: | N/A | 34 | 34 | 34 | N/A |
| Maximum treadmill height | Treadmill/stability bases: | N/A | 11 | 11 | 11 | N/A |
| Maximum working load (lb) | | 50 | 75 | 150 | 200 | 250 |

*This measurement is taken at the lowest frame and MPS height, and the highest frame and MPS height.

Pelvic support

| User dimensions (inches) Indicated by the red arrow below | | | |
|---|-----------|------------|----|
| Key user dimension: half of thigh circumference measured at the groin. | | | |
| K556 small | K557 med. | K558 large | |
| 4–8 | 6–10 | 8–14 | |
|  | | | |
| Product dimensions (inches) | | | |
| K556 small | K557 med. | K558 large | |
| a | 6½ | 8½ | 12 |
| b | 12 | 16 | 20 |
| c | 7 | 9 | 12 |

Hip positioner

| Dimensions (inches) – without pad | | |
|---|------------|------------|
|  | | |
| | K513 small | K533 large |
| a Crotch length | 2½ | 3½ |
| b Width | 1¼ | 1¾ |
| c Length | 9 | 11 |



All five sizes are available in five colors: red, blue, gray, lime and pink.