


# E-Pacer 2018

1. **Choose** your color.
2. **Select** additional accessories as your needs require.
3. **Add** up the total cost and place your order.
4. **Order** online in minutes at **Rifton.com** or see the back of this catalog for other ways to order.
5. **Call** us any time at **800.571.8198** if you need assistance. We're always here to help.




		K660		
Required components	<b>Frame</b> includes one battery, one battery charger and charger power cord (international power cords available on request)	K661	\$4,085 <input type="checkbox"/>	
	<b>Color</b>	CHOOSE	Red <input type="checkbox"/>	
			Blue <input type="checkbox"/>	
			Gray <input type="checkbox"/>	
			Lime <input type="checkbox"/>	
Pink <input type="checkbox"/>				
Optional accessories	<b>Scale</b> (must be chosen at the time of the original purchase. It cannot be added later.)	K662	\$965 <input type="checkbox"/>	
	<b>Odometer</b>	K663	\$160 <input type="checkbox"/>	
	<b>Front handle</b>	K665	\$310 <input type="checkbox"/>	
	<b>Switch pole</b>	K666	\$205 <input type="checkbox"/>	
	<b>Arm supports</b> (pair)	CHOOSE	Small arm prompts	K305 \$450 <input type="checkbox"/>
			Large arm prompts	K315 \$520 <input type="checkbox"/>
			Arm platforms	K397 \$335 <input type="checkbox"/>
			Arm platforms with handgrips	K396 \$385 <input type="checkbox"/>
	<b>Hip positioner</b>	CHOOSE	Small	K667 \$280 <input type="checkbox"/>
			Large	K668 \$310 <input type="checkbox"/>
	<b>Hip positioner pad</b>	CHOOSE	Small	K541 \$98 <input type="checkbox"/>
			Large	K542 \$130 <input type="checkbox"/>
	<b>Pelvic support</b>	CHOOSE	Small	K306 \$230 <input type="checkbox"/>
			Medium	K316 \$255 <input type="checkbox"/>
			Large	K326 \$270 <input type="checkbox"/>
	<b>Thigh straps</b>	CHOOSE	Narrow	K308 \$68 <input type="checkbox"/>
			Wide	K338 \$68 <input type="checkbox"/>
	<b>Ankle prompts</b>		K684	\$190 <input type="checkbox"/>
	<b>Additional battery</b>		K307	\$125 <input type="checkbox"/>
	<b>Additional thigh straps</b>	CHOOSE	Narrow single	K318 \$34 <input type="checkbox"/>
Narrow pair			K308 \$68 <input type="checkbox"/>	
Wide single			K328 \$34 <input type="checkbox"/>	
Wide pair			K338 \$68 <input type="checkbox"/>	
<b>Removable back belt</b>		K327	\$53 <input type="checkbox"/>	
<b>Rifton accessories tote</b>		K522	\$88 <input type="checkbox"/>	
<b>Total</b>		<b>K660</b>	<b>\$_____</b>	

# Dimensions of the E-Pacer

User dimensions (inches)	K660
Girth	22 min–60 max
 <p><b>Key user dimension: girth</b></p> <p><b>Important:</b> User's weight or height must not exceed: 350 lbs or 77 inches</p>	
Item dimensions (inches)	K660 E-Pacer
Overall length	44¾
Overall width	31½–51¾
Overall height	32¾–59
E-Pacer frame weight (lbs)	72
Turning diameter	50
Floor to top of base leg for furniture clearance	8¾
Min. user armpit height (when standing or ambulating)	31½
Max. user height (when standing or ambulating)	77
Max. working load (lbs)	350

## Pelvic support

User dimensions (inches) Indicated by the red arrow below			
Key user dimension: half of thigh circumference.			
	K306 small	K316 med.	K326 large
	6–10	8–14	10–18
			
Product dimensions (inches)			
	K306 small	K316 med.	K326 large
<b>a</b>	8½	12	13½
<b>b</b>	16	20	25
<b>c</b>	9	12	14

The E-Pacer frame comes in your choice of five colors.



Red



Blue



Gray



Lime



Pink

## Hip positioner

Dimensions (inches) – without pad		
	K667 small	K668 large
<b>a</b> Crotch length	2½	3½
<b>b</b> Width	1¼	1¾
<b>c</b> Length	9	11