

2021 order form


Use dimension chart to select appropriate size.




| | | K660 | |
|--|--------|------------------------------|-------------------------------------|
| * Frame includes one battery, one battery charger and charger power cord. (International power cords available on request.) | | K661 | \$4,460 <input type="checkbox"/> |
| * Color | | CHOOSE | Red <input type="checkbox"/> |
| | | | Blue <input type="checkbox"/> |
| | | | Gray <input type="checkbox"/> |
| | | | Lime <input type="checkbox"/> |
| | | | Pink <input type="checkbox"/> |
| Scale (must be chosen at the time of the original purchase. It cannot be added later.) | | K662 | \$1,060 <input type="checkbox"/> |
| Odometer | | K663 | \$175 <input type="checkbox"/> |
| Front handle | | K665 | \$335 <input type="checkbox"/> |
| Switch pole | | K666 | \$225 <input type="checkbox"/> |
| Arm supports (pair) | CHOOSE | Small arm prompts | K305 \$475 <input type="checkbox"/> |
| | | Large arm prompts | K315 \$545 <input type="checkbox"/> |
| | | Arm platforms | K397 \$350 <input type="checkbox"/> |
| | | Arm platforms with handgrips | K396 \$400 <input type="checkbox"/> |
| Hip positioner | CHOOSE | Small | K667 \$295 <input type="checkbox"/> |
| | | Large | K668 \$325 <input type="checkbox"/> |
| Hip positioner pad | CHOOSE | Small | K541 \$100 <input type="checkbox"/> |
| | | Large | K542 \$135 <input type="checkbox"/> |
| Pelvic support | CHOOSE | Small | K306 \$240 <input type="checkbox"/> |
| | | Medium | K316 \$270 <input type="checkbox"/> |
| | | Large | K326 \$285 <input type="checkbox"/> |
| Thigh straps | CHOOSE | Narrow | K308 \$100 <input type="checkbox"/> |
| | | Wide | K338 \$100 <input type="checkbox"/> |
| Ankle prompts | | K684 | \$200 <input type="checkbox"/> |
| Additional battery | | K307 | \$135 <input type="checkbox"/> |
| Additional thigh straps | CHOOSE | Narrow single | K318 \$50 <input type="checkbox"/> |
| | | Narrow pair | K308 \$100 <input type="checkbox"/> |
| | | Wide single | K328 \$50 <input type="checkbox"/> |
| | | Wide pair | K338 \$100 <input type="checkbox"/> |
| Removable back belt | | K327 | \$60 <input type="checkbox"/> |
| Rifton accessories tote | | K522 | \$120 <input type="checkbox"/> |
| Total | | K660 | \$ _____ |

* Required


E-Pacer dimensions

| User dimensions (inches) | | K660 |
|---|----------------------------------|--------------|
| Girth | 22 min–60 max | |
|  | Key user dimension: girth | |
| Important: User's weight or height must not exceed: 350 lb or 77 inches | | |
| Item dimensions (inches) | | K660 E-Pacer |
| Overall length | 44¾ | |
| Overall width | 31½–51¾ | |
| Overall height | 32¾–59 | |
| E-Pacer frame weight (lb) | 72 | |
| Turning diameter | 50 | |
| Floor to top of base leg for furniture clearance | 8¾ | |
| Min. user armpit height (when standing or ambulating) | 31½ | |
| Max. user height (when standing or ambulating) | 77 | |
| Max. working load (lb) | 350 | |

Pelvic support

| User dimensions (inches) Indicated by the red arrow below | | | |
|---|-----------|------------|-----|
| Key user dimension: half of thigh circumference measured at the groin. | | | |
| K306 small | K316 med. | K326 large | |
| 6–10 | 8–14 | 10–18 | |
|  | | | |
| Product dimensions (inches) | | | |
| K306 small | K316 med. | K326 large | |
| a | 8½ | 12 | 13½ |
| b | 16 | 20 | 25 |
| c | 9 | 12 | 14 |


Hip positioner

| Product dimensions (inches)–without pad | | |
|--|------------|------------|
|  | | |
| | K667 small | K668 large |
| a Crotch length | 2½ | 3½ |
| b Width | 1¼ | 1¾ |
| c Length | 9 | 11 |

The E-Pacer comes in your choice of five colors.



Red



Blue



Gray



Lime



Pink