

Mobile Stander 2018

- 1. Measure** your client to determine the size you need. (See dimension chart for help.)
- 2. Select** the appropriate size and color.
- 3. Choose** the accessories you need.
- 4. Add** up the total cost.
- 5. Order** online in minutes at **Rifton.com** or see the back of this catalog for other ways to order.
- 6. Call us** any time at **800.571.8198** if you need assistance. We're always here to help.



The **conversion kit** converts our smallest stander, the K110, to the next size up, the K130, for less than one quarter the cost of a new stander. Request from customer service.

Conversion kit K116 \$650



Mobile Stander frame includes body support with padding, seat pad, straps and fleece.

Specify color: pink, blue or lime

	K110 mini		K130 small		K150 medium		K170 large	
CHOOSE	K101	\$1,760 <input type="checkbox"/>	K103	\$1,960 <input type="checkbox"/>	K105	\$2,390 <input type="checkbox"/>	K107	\$2,920 <input type="checkbox"/>
		PNK <input type="checkbox"/>		PNK <input type="checkbox"/>		PNK <input type="checkbox"/>		PNK <input type="checkbox"/>
		BLU <input type="checkbox"/>		BLU <input type="checkbox"/>		BLU <input type="checkbox"/>		BLU <input type="checkbox"/>
		LIM <input type="checkbox"/>		LIM <input type="checkbox"/>		LIM <input type="checkbox"/>		LIM <input type="checkbox"/>

		K110 mini		K130 small		K150 medium		K170 large			
		Item	Price	Item	Price	Item	Price	Item	Price		
Optional accessories	Wheels (pair) (includes brakes)	K113	\$1,045 <input type="checkbox"/>	K113	\$1,045 <input type="checkbox"/>	K115	\$1,185 <input type="checkbox"/>	K117	\$1,425 <input type="checkbox"/>		
	Mounting bracket	K133	\$57 <input type="checkbox"/>	K133	\$57 <input type="checkbox"/>	K173	\$67 <input type="checkbox"/>	K173	\$67 <input type="checkbox"/>		
	Arm prompts (pair) (require mounting bracket)	K512	\$455 <input type="checkbox"/>	K512	\$455 <input type="checkbox"/>	K512	\$455 <input type="checkbox"/>	K512	\$455 <input type="checkbox"/>		
	Hand loops (require mounting bracket)	CHOOSE	Single	K578	\$90 <input type="checkbox"/>	K578	\$90 <input type="checkbox"/>	K578	\$90 <input type="checkbox"/>	K578	\$90 <input type="checkbox"/>
			Pair	K579	\$180 <input type="checkbox"/>	K579	\$180 <input type="checkbox"/>	K579	\$180 <input type="checkbox"/>	K579	\$180 <input type="checkbox"/>
	Communication tray (requires mounting bracket)	K518	\$345 <input type="checkbox"/>	K518	\$345 <input type="checkbox"/>	K518	\$345 <input type="checkbox"/>	K518	\$345 <input type="checkbox"/>		
	Additional seat pad	CHOOSE	Small	N/A	N/A	K135	\$130 <input type="checkbox"/>	K135	\$130 <input type="checkbox"/>		
			Large	N/A	N/A	K155	\$140 <input type="checkbox"/>	K155	\$140 <input type="checkbox"/>		
	Sandals with wedges (pair)	CHOOSE	Small	K114	\$98 <input type="checkbox"/>	K114	\$98 <input type="checkbox"/>	N/A	N/A		
			Medium	N/A	N/A	K154	\$105 <input type="checkbox"/>	N/A			
Large			N/A	N/A	N/A	N/A	K174	\$110 <input type="checkbox"/>			
Total		K110	\$ _____	K130	\$ _____	K150	\$ _____	K170	\$ _____		

User dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Height	25-32	30-40	37-52	48-65
Key user dimension: height Select the appropriate Mobile Stander by the user's overall height. Choose the model that allows for growth. Important: User's weight must not exceed the maximum working load.				
Item dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Overall length and width	25x25	25x25½	32x29	39x33½
Overall height	22½-28	23½-30	31-39	37½-49
Width without big wheels	19½	21	24	28½
Wheel diameter	20	20	27	35
Width of body support	7-10	7-10	9-12	10-14
Height of body support	16½-22	22-29	29-38	36-48
Item weight (lbs)	28	30	40	65
Max. working load (lbs)	50	50	120	175