


2019 order form

Use dimension chart to select appropriate size.

User dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Height	25–32	30–40	37–52	48–65
 <p>Key user dimension: height Select the appropriate Mobile Stander by the user's overall height. Choose the model that allows for growth.</p> <p>Important: User's weight must not exceed the maximum working load.</p>				
	Item dimensions (inches)	K110 mini	K130 small	K150 med.
Overall length and width	25x25	25x25½	32x29	39x33½
Overall height	22½–28	23½–30	31–39	37½–49
Width without large wheels	19½	19½	24½	28¾
Wheel diameter	20	20	27	35
Width of body support	7–10	7–10	9–12	10–14
Height of body support	16½–22	22–29	29–38	36–48
Item weight (lbs)	28	30	40	65
Max. working load (lbs)	50	50	120	175



The **conversion kit** converts our smallest stander, the K110, to the next size up, the K130, for less than one quarter the cost of a new stander. Request from customer service.

Conversion kit K116 \$675



* **Mobile Stander** frame includes body support with padding, seat pad, straps and fleece.

Specify color: pink, blue or lime

		K110 mini	K130 small	K150 medium	K170 large	
CHOOSE	K101	\$1,830 <input type="checkbox"/>	K103 \$2,040 <input type="checkbox"/>	K105 \$2,485 <input type="checkbox"/>	K107 \$3,035 <input type="checkbox"/>	
	Pink	<input type="checkbox"/>	Pink <input type="checkbox"/>	Pink <input type="checkbox"/>	Pink <input type="checkbox"/>	
	Blue	<input type="checkbox"/>	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>	
	Lime	<input type="checkbox"/>	Lime <input type="checkbox"/>	Lime <input type="checkbox"/>	Lime <input type="checkbox"/>	
Wheels (pair) (includes brakes)		K113 \$1,085 <input type="checkbox"/>	K113 \$1,085 <input type="checkbox"/>	K115 \$1,230 <input type="checkbox"/>	K117 \$1,480 <input type="checkbox"/>	
Mounting bracket		K133 \$59 <input type="checkbox"/>	K133 \$59 <input type="checkbox"/>	K173 \$70 <input type="checkbox"/>	K173 \$70 <input type="checkbox"/>	
Arm prompts (pair) (require mounting bracket)		K512 \$475 <input type="checkbox"/>	K512 \$475 <input type="checkbox"/>	K512 \$475 <input type="checkbox"/>	K512 \$475 <input type="checkbox"/>	
Hand loops (require mounting bracket)	CHOOSE	Single	K578 \$94 <input type="checkbox"/>	K578 \$94 <input type="checkbox"/>	K578 \$94 <input type="checkbox"/>	K578 \$94 <input type="checkbox"/>
	Pair	K579 \$185 <input type="checkbox"/>	K579 \$185 <input type="checkbox"/>	K579 \$185 <input type="checkbox"/>	K579 \$185 <input type="checkbox"/>	
Communication tray (requires mounting bracket)		K518 \$360 <input type="checkbox"/>	K518 \$360 <input type="checkbox"/>	K518 \$360 <input type="checkbox"/>	K518 \$360 <input type="checkbox"/>	
Additional seat pad	CHOOSE	Small	N/A	N/A	K135 \$135 <input type="checkbox"/>	K135 \$135 <input type="checkbox"/>
	Large	N/A	N/A	K155 \$145 <input type="checkbox"/>	K155 \$145 <input type="checkbox"/>	
Sandals with wedges (pair)	CHOOSE	Small	K114 \$100 <input type="checkbox"/>	K114 \$100 <input type="checkbox"/>	N/A	N/A
	Medium	N/A	N/A	K154 \$110 <input type="checkbox"/>	N/A	
	Large	N/A	N/A	N/A	K174 \$115 <input type="checkbox"/>	
Total		K110 \$ _____	K130 \$ _____	K150 \$ _____	K170 \$ _____	

* Required