

2020 order form

Use dimension chart to select appropriate size.



The **conversion kit** converts the K110 mini stander to the K130 small. It includes the deck and complete body support of the K130. Request from customer service.

Conversion kit K116 \$710




* **Mobile Stander** frame includes body support with padding, seat pad, straps and fleece.

Specify color:
pink, blue or lime

		K110 mini		K130 small		K150 medium		K170 large	
CHOOSE		K101	\$1,900 <input type="checkbox"/>	K103	\$2,150 <input type="checkbox"/>	K105	\$2,600 <input type="checkbox"/>	K107	\$3,100 <input type="checkbox"/>
		Pink	<input type="checkbox"/>	Pink	<input type="checkbox"/>	Pink	<input type="checkbox"/>	Pink	<input type="checkbox"/>
		Blue	<input type="checkbox"/>	Blue	<input type="checkbox"/>	Blue	<input type="checkbox"/>	Blue	<input type="checkbox"/>
		Lime	<input type="checkbox"/>	Lime	<input type="checkbox"/>	Lime	<input type="checkbox"/>	Lime	<input type="checkbox"/>
Wheels (pair) (includes brakes)		K113	\$1,140 <input type="checkbox"/>	K113	\$1,140 <input type="checkbox"/>	K115	\$1,290 <input type="checkbox"/>	K117	\$1,550 <input type="checkbox"/>
Mounting bracket		K133	\$65 <input type="checkbox"/>	K133	\$65 <input type="checkbox"/>	K173	\$75 <input type="checkbox"/>	K173	\$75 <input type="checkbox"/>
Arm prompts (pair) (require mounting bracket)		K512	\$485 <input type="checkbox"/>	K512	\$485 <input type="checkbox"/>	K512	\$485 <input type="checkbox"/>	K512	\$485 <input type="checkbox"/>
Hand loops (require mounting bracket)	CHOOSE Single	K578	\$80 <input type="checkbox"/>	K578	\$80 <input type="checkbox"/>	K578	\$80 <input type="checkbox"/>	K578	\$80 <input type="checkbox"/>
	Pair	K579	\$190 <input type="checkbox"/>	K579	\$190 <input type="checkbox"/>	K579	\$190 <input type="checkbox"/>	K579	\$190 <input type="checkbox"/>
Communication tray (requires mounting bracket)		K518	\$370 <input type="checkbox"/>	K518	\$370 <input type="checkbox"/>	K518	\$370 <input type="checkbox"/>	K518	\$370 <input type="checkbox"/>
Additional seat pad	CHOOSE Small	N/A		N/A		K135	\$140 <input type="checkbox"/>	K135	\$140 <input type="checkbox"/>
	Large	N/A		N/A		K155	\$150 <input type="checkbox"/>	K155	\$150 <input type="checkbox"/>
Sandals with wedges (pair)	CHOOSE Small	K114	\$105 <input type="checkbox"/>	K114	\$105 <input type="checkbox"/>	N/A		N/A	
	Medium	N/A		N/A		K154	\$115 <input type="checkbox"/>	N/A	
	Large	N/A		N/A		N/A		K174	\$120 <input type="checkbox"/>
Total		K110	\$ _____	K130	\$ _____	K150	\$ _____	K170	\$ _____

* Required

User dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Height	25–32	30–40	37–52	48–65
 Key user dimension: height Select the appropriate Mobile Stander by the user's overall height. Choose the model that allows for growth.	Important: User's weight must not exceed the maximum working load.			
Item dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Overall length and width	25x25	25x25½	32x29	39x33½
Overall height	22½–28	23½–30	31–39	37½–49
Width without large wheels	19½	19½	24½	28¾
Wheel diameter	20	20	27	35
Width of body support	7–10	7–10	9–12	10–14
Height of body support	16½–22	22–29	29–38	36–48
Item weight (lbs)	28	30	40	65
Max. working load (lbs)	50	50	120	175