

2021 order form

Use dimension chart to select appropriate size.



The **conversion kit** converts the K110 mini stander to the K130 small. It includes the deck and complete body support of the K130. Request from customer service.

Conversion kit K116 \$710




* **Mobile Stander** frame includes body support with padding, seat pad, straps and fleece.

Specify color:
pink, blue or lime

		K110 mini	K130 small	K150 medium	K170 large			
CHOOSE	K101	\$1,900 <input type="checkbox"/>	K103 \$2,150 <input type="checkbox"/>	K105 \$2,600 <input type="checkbox"/>	K107 \$3,100 <input type="checkbox"/>			
	Pink	<input type="checkbox"/>	Pink <input type="checkbox"/>	Pink <input type="checkbox"/>	Pink <input type="checkbox"/>			
	Blue	<input type="checkbox"/>	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>			
	Lime	<input type="checkbox"/>	Lime <input type="checkbox"/>	Lime <input type="checkbox"/>	Lime <input type="checkbox"/>			
Wheels (pair) (include brakes)	K113	\$1,140 <input type="checkbox"/>	K113 \$1,140 <input type="checkbox"/>	K115 \$1,290 <input type="checkbox"/>	K117 \$1,550 <input type="checkbox"/>			
Mounting bracket	K133	\$65 <input type="checkbox"/>	K133 \$65 <input type="checkbox"/>	K173 \$75 <input type="checkbox"/>	K173 \$75 <input type="checkbox"/>			
Arm prompts (pair) (require mounting bracket)	K512	\$485 <input type="checkbox"/>	K512 \$485 <input type="checkbox"/>	K512 \$485 <input type="checkbox"/>	K512 \$485 <input type="checkbox"/>			
Hand loops (require mounting bracket)	CHOOSE	Single	K578 \$80 <input type="checkbox"/>	K578 \$80 <input type="checkbox"/>	K578 \$80 <input type="checkbox"/>	K578 \$80 <input type="checkbox"/>		
		Pair	K579 \$190 <input type="checkbox"/>	K579 \$190 <input type="checkbox"/>	K579 \$190 <input type="checkbox"/>	K579 \$190 <input type="checkbox"/>		
Communication tray (requires mounting bracket)	K518	\$370 <input type="checkbox"/>	K518 \$370 <input type="checkbox"/>	K518 \$370 <input type="checkbox"/>	K518 \$370 <input type="checkbox"/>			
Additional seat pad	CHOOSE	Small	N/A	N/A	K135 \$140 <input type="checkbox"/>	K135 \$140 <input type="checkbox"/>		
		Large	N/A	N/A	K155 \$150 <input type="checkbox"/>	K155 \$150 <input type="checkbox"/>		
Sandals with wedges (pair)	CHOOSE	Small	K114 \$105 <input type="checkbox"/>	K114 \$105 <input type="checkbox"/>	N/A	N/A		
		Medium	N/A	N/A	K154 \$115 <input type="checkbox"/>	N/A		
		Large	N/A	N/A	N/A	K174 \$120 <input type="checkbox"/>		
Total	K110	\$ _____	K130	\$ _____	K150	\$ _____	K170	\$ _____

* Required

User dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Height	25–32	30–40	37–52	48–65
 Key user dimension: height Select the appropriate Mobile Stander by the user's overall height. Choose the model that allows for growth.	Important: User's weight must not exceed the maximum working load.			
Item dimensions (inches)	K110 mini	K130 small	K150 med.	K170 large
Overall length and width	25x25	25x25½	32x29	39x33½
Overall height	22½–28	23½–30	31–39	37½–49
Width without large wheels	19½	19½	24½	28¾
Wheel diameter	20	20	27	35
Width of body support	7–10	7–10	9–12	10–14
Height of body support	16½–22	22–29	29–38	36–48
Item weight (lb)	28	30	40	65
Max. working load (lb)	50	50	120	175