



Pacer Gait Trainers 2018

- 1. Measure** your client to determine the size you need.
(See dimension chart for help.)
- 2. Select** frame and color.
- 3. Choose** the accessories you need.
- 4. Add** up the total cost.
- 5. Order** online in minutes at **Rifton.com** or see the back of this catalog for other ways to order.
- 6. Call** us any time at **800.571.8198** if you need assistance.
We're always here to help.

User dimensions (inches)	K509 mini	K501 small
Elbow height	15½ – 20½	18½ – 27½
 <p>Key user dimension: elbow height Measure the vertical distance from bent elbow to floor while the user is standing upright. Choose the gait trainer that allows for growth.</p> <p>Important: Make sure that chest prompt width allows for growth. User's weight must not exceed the maximum working load.</p>		
Item dimensions (inches)	K509 mini	K501 small
Arm prompt height	15½ – 20½	18½ – 27½
Overall width	20½	22½
Overall length	22½	27
Chest prompt height (top edge)	17½ – 22½	22 – 30½
Frame height	12½	15½ – 19½
Frame weight (lbs)	7½	11½
Max. working load (lbs)	50	75



	K509 mini		K501 small	
Frame	K510	\$680 <input type="checkbox"/>	K511	\$920 <input type="checkbox"/>
Specify color: Red, blue, lime or pink	RED <input type="checkbox"/>	PNK <input type="checkbox"/>	RED <input type="checkbox"/>	PNK <input type="checkbox"/>
	BLU <input type="checkbox"/>	LIM <input type="checkbox"/>	BLU <input type="checkbox"/>	LIM <input type="checkbox"/>

		K509 mini		K501 small		
Optional accessories	Arm prompts (pair)	K512	\$455 <input type="checkbox"/>	K512	\$455 <input type="checkbox"/>	
	Hand loops	CHOOSE	Single	K578 \$90 <input type="checkbox"/>	K578 \$90 <input type="checkbox"/>	
			Pair	K579 \$180 <input type="checkbox"/>	K579 \$180 <input type="checkbox"/>	
	Chest prompt	CHOOSE	Small	K514 \$510 <input type="checkbox"/>	K514 \$510 <input type="checkbox"/>	
			Medium	N/A	K524 \$550 <input type="checkbox"/>	
	Pelvic support (includes handholds)	CHOOSE	Small	K556 \$385 <input type="checkbox"/>	K556 \$385 <input type="checkbox"/>	
			Medium	N/A	K557 \$405 <input type="checkbox"/>	
	Hip positioner (includes handholds)		K513	\$385 <input type="checkbox"/>	K513	\$385 <input type="checkbox"/>
	Hip positioner pad		K541	\$98 <input type="checkbox"/>	K541	\$98 <input type="checkbox"/>
	Thigh prompts (pair)		N/A		K515	\$380 <input type="checkbox"/>
	Ankle prompts (pair)		K529	\$150 <input type="checkbox"/>	K516	\$160 <input type="checkbox"/>
	Additional handholds (pair)		K519	\$105 <input type="checkbox"/>	K519	\$105 <input type="checkbox"/>
	Communication tray		K518	\$345 <input type="checkbox"/>	K518	\$345 <input type="checkbox"/>
Attendant guide bar		K589	\$130 <input type="checkbox"/>	K589	\$130 <input type="checkbox"/>	
Accessories tote		K522	\$88 <input type="checkbox"/>	K522	\$88 <input type="checkbox"/>	
Total		K509	\$ _____	K501	\$ _____	