


2019 order form

Use dimension chart to select appropriate size.



User dimensions (inches)	E930 small	E940 med.	E950 large
Height	25 – 48	44 – 59	57 – 72
 <p>Key user dimension: height Select the appropriate Prone Stander by the user's overall height. Choose the model that allows for growth.</p> <p>Important: User's weight must not exceed the maximum working load.</p>			
	Item dimensions (inches)	E930 small	E940 med.
Board length	23½ (16 without kneeboard)	30	39
Board width (top / bottom)	12½ / 12½	11 / 16	14 / 19
Top of board to footboard	20 – 35	33 – 43½	44 – 57
Height when horizontal	29½	29½	31½
Height when vertical	38½	45½	61
Distance between laterals	8 – 16	8 – 16	12 – 22
Base width	23	26	30
Base length	35	39½	44
Tray inside width (widest point)	22½	22½	29½
Tray inside length	17½	17½	17½
Max. working load (lbs)	100	150	200

* **Prone Stander** includes a pair of trunk laterals, hip strap, padding and footboard.

Specify color: red or blue

		E930 small	E940 medium	E950 large		
CHOOSE	E931	\$1,990 <input type="checkbox"/>	E941 \$2,305 <input type="checkbox"/>	E951 \$2,820 <input type="checkbox"/>		
	Red	<input type="checkbox"/>	Red <input type="checkbox"/>	Red <input type="checkbox"/>		
	Blue	<input type="checkbox"/>	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>		
Additional trunk laterals (pair)	N/A	E614 \$200 <input type="checkbox"/>	E614 \$200 <input type="checkbox"/>			
Knee laterals (pair)	E942 \$185 <input type="checkbox"/>	E942 \$185 <input type="checkbox"/>	E942 \$185 <input type="checkbox"/>			
Additional hip strap	N/A	E946 \$70 <input type="checkbox"/>	E956 \$80 <input type="checkbox"/>			
Tray (see dimension chart for size)	E947 \$555 <input type="checkbox"/>	E947 \$555 <input type="checkbox"/>	E957 \$640 <input type="checkbox"/>			
Round abductor	E933 \$95 <input type="checkbox"/>	E933 \$95 <input type="checkbox"/>	E933 \$95 <input type="checkbox"/>			
Collar for round abductor	E934 \$75 <input type="checkbox"/>	E934 \$75 <input type="checkbox"/>	E934 \$75 <input type="checkbox"/>			
Round abductor with hip stabilizer	E655 \$270 <input type="checkbox"/>	E655 \$270 <input type="checkbox"/>	E655 \$270 <input type="checkbox"/>			
Adjustable abduction wedge	N/A	E945 \$225 <input type="checkbox"/>	E945 \$225 <input type="checkbox"/>			
Adjustable abduction wedge with hip stabilizer	N/A	E605 \$285 <input type="checkbox"/>	E605 \$285 <input type="checkbox"/>			
Sandals (pair)	CHOOSE	Small	E422 \$85 <input type="checkbox"/>	E422 \$85 <input type="checkbox"/>	N/A	
		Medium	E433 \$94 <input type="checkbox"/>	E433 \$94 <input type="checkbox"/>	E433 \$94 <input type="checkbox"/>	
		Large	N/A	N/A	E432 \$97 <input type="checkbox"/>	
Wedges (pair, requires sandals)	CHOOSE	Small	E424 \$61 <input type="checkbox"/>	E424 \$61 <input type="checkbox"/>	N/A	
		Medium	E435 \$67 <input type="checkbox"/>	E435 \$67 <input type="checkbox"/>	E435 \$67 <input type="checkbox"/>	
		Large	N/A	N/A	E434 \$70 <input type="checkbox"/>	
Hand Anchor (For more hand/arm support options see pp. 98 – 101.)	K821 \$145 <input type="checkbox"/>	K821 \$145 <input type="checkbox"/>	K821 \$145 <input type="checkbox"/>			
Total	E930	\$ _____	E940	\$ _____	E950	\$ _____

* Required