

# Large Activity Chair 2018



			Standard R860 large		Hi/lo R870 large		
Required components	<b>Seat and back</b> includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE Adjustable	R861	\$605 <input type="checkbox"/>	R871	\$1,125 <input type="checkbox"/>	
		Adjustable w/ back spring	R862	\$700 <input type="checkbox"/>	R872	\$1,215 <input type="checkbox"/>	
	<b>Base</b>	CHOOSE Without spring	R873	\$370 <input type="checkbox"/>	R883	\$2,500 <input type="checkbox"/>	
		With spring	R874	\$440 <input type="checkbox"/>	N/A		
	<b>Arm supports</b> (pair)	CHOOSE Armrests	R868	\$215 <input type="checkbox"/>	R868	\$215 <input type="checkbox"/>	
		Forearm prompts	R879	\$560 <input type="checkbox"/>	R879	\$560 <input type="checkbox"/>	
	<b>Adjustable legs</b>	CHOOSE	Short	R808	\$59 <input type="checkbox"/>	N/A	
			Long	R894	\$70 <input type="checkbox"/>	N/A	
			Short with casters	R895	\$260 <input type="checkbox"/>	N/A	
			Long with casters	R896	\$285 <input type="checkbox"/>	N/A	
	<b>Pads</b>	CHOOSE	Tan	R875	\$225 <input type="checkbox"/>	R875	\$225 <input type="checkbox"/>
			Blue	R876	\$225 <input type="checkbox"/>	R876	\$225 <input type="checkbox"/>
			Red	R877	\$225 <input type="checkbox"/>	R877	\$225 <input type="checkbox"/>
			Pink	R977	\$225 <input type="checkbox"/>	R977	\$225 <input type="checkbox"/>
			Purple	R978	\$225 <input type="checkbox"/>	R978	\$225 <input type="checkbox"/>
			Green	R979	\$225 <input type="checkbox"/>	R979	\$225 <input type="checkbox"/>
Optional accessories	<b>Headrest</b>	CHOOSE Contoured	R897	\$210 <input type="checkbox"/>	R897	\$210 <input type="checkbox"/>	
		Flat	R898	\$195 <input type="checkbox"/>	R898	\$195 <input type="checkbox"/>	
		Adjustable winged	R881	\$260 <input type="checkbox"/>	R881	\$260 <input type="checkbox"/>	
	<b>Laterals</b> (pair)	CHOOSE	Small	R828	\$175 <input type="checkbox"/>	R828	\$175 <input type="checkbox"/>
			Large	R838	\$190 <input type="checkbox"/>	R838	\$190 <input type="checkbox"/>
			Large with strap	R839	\$225 <input type="checkbox"/>	R839	\$225 <input type="checkbox"/>
	<b>Push handles</b> (recommended for all chairs with wheels)		R878	\$135 <input type="checkbox"/>	Included		
	<b>Tray</b> (requires armrests)		R669	\$250 <input type="checkbox"/>	R669	\$250 <input type="checkbox"/>	
	<b>Handhold</b>		R649	\$21 <input type="checkbox"/>	R649	\$21 <input type="checkbox"/>	
	<b>Chest strap</b> (wide)		R807	\$82 <input type="checkbox"/>	R807	\$82 <input type="checkbox"/>	
	<b>Butterfly harness</b>	CHOOSE	Regular	R889	\$130 <input type="checkbox"/>	R889	\$130 <input type="checkbox"/>
			Slim cut	R826	\$130 <input type="checkbox"/>	R826	\$130 <input type="checkbox"/>
	<b>Thigh belt</b>		R813	\$82 <input type="checkbox"/>	R813	\$82 <input type="checkbox"/>	
	<b>Pelvic harness</b>		R846	\$155 <input type="checkbox"/>	R846	\$155 <input type="checkbox"/>	
	<b>Hip guides</b> (pair)		R867	\$185 <input type="checkbox"/>	R867	\$185 <input type="checkbox"/>	
	<b>Abductor</b>	CHOOSE	Small	R884	\$110 <input type="checkbox"/>	R884	\$110 <input type="checkbox"/>
			Large	R891	\$130 <input type="checkbox"/>	R891	\$130 <input type="checkbox"/>
	<b>Adductors</b> (pair)		R892	\$175 <input type="checkbox"/>	R892	\$175 <input type="checkbox"/>	
	<b>Leg prompt</b>		R859	\$205 <input type="checkbox"/>	R859	\$205 <input type="checkbox"/>	
	<b>Footboard</b> (recommended for chairs with long legs)		R857	\$410 <input type="checkbox"/>	Included		
	<b>Footboard lift</b>		R818	\$210 <input type="checkbox"/>	R818	\$210 <input type="checkbox"/>	
	<b>Ankle straps</b> (pair)		R855	\$83 <input type="checkbox"/>	R855	\$83 <input type="checkbox"/>	
<b>Sandals</b> (pair)		R964	\$94 <input type="checkbox"/>	R964	\$94 <input type="checkbox"/>		
<b>Wedges</b> (pair, requires sandals)		R699	\$64 <input type="checkbox"/>	R699	\$64 <input type="checkbox"/>		
<b>Lumbar and seat support kit</b>		R888	\$77 <input type="checkbox"/>	R888	\$77 <input type="checkbox"/>		
<b>Backrest filler pad</b>		R856	\$110 <input type="checkbox"/>	R856	\$110 <input type="checkbox"/>		
<b>Whitmyer headrest adapter</b>		R804	\$40 <input type="checkbox"/>	R804	\$40 <input type="checkbox"/>		
<b>Total</b>		<b>R860</b>	<b>\$ _____</b>	<b>R870</b>	<b>\$ _____</b>		

# Room to grow

## Large chair at its smallest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 9½"
- 3 Seat depth 15"
- 4 Seat width 11"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 15" (with footboard lift 11")



## Large chair at its largest dimensions


- 1 Backrest height 24"
- 2 Distance between laterals 14"
- 3 Seat depth 20"
- 4 Seat width 17"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 20"




**Large standard base**  
R860 Rifton Activity Chair



**Large hi/lo base**  
R870 Rifton Activity Chair

User dimensions (inches)	
Height	50–74
 <p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	
	R860 large Standard base
Frame width	short legs: 26¾ long legs: 26¾
	short legs w/ casters: 27¾ long legs w/ casters: 27¾
Seat height above floor	short legs: 16½–20½ long legs: 20–24
	short legs w/ casters: 16½–20½ long legs w/ casters: 20–24
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Seat height above footboard with footboard lift	15–20 11–16
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lbs)	250

User dimensions (inches)	
Height	50–74
 <p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	
	R870 large Hi/lo base
Frame width	29¾
Seat height above floor	16–24
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Seat height above footboard with footboard lift	15–20 11–16
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lbs)	250