

2019 order form for large

Use dimension chart to select appropriate size.



		Standard R860 large		Hi/lo R870 large		
* Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R861	\$630 <input type="checkbox"/>	R871	\$1,200 <input type="checkbox"/>
		Adjustable w/ back spring	R862	\$730 <input type="checkbox"/>	R872	\$1,300 <input type="checkbox"/>
* Base	CHOOSE	Without spring	R873	\$385 <input type="checkbox"/>	R883	\$2,600 <input type="checkbox"/>
		With spring	R874	\$460 <input type="checkbox"/>	N/A	
* Arm supports (pair)	CHOOSE	Armrests	R868	\$225 <input type="checkbox"/>	R868	\$225 <input type="checkbox"/>
		Forearm prompts	R879	\$580 <input type="checkbox"/>	R879	\$580 <input type="checkbox"/>
* Adjustable legs	CHOOSE	Short	R808	\$61 <input type="checkbox"/>	N/A	
		Long	R894	\$73 <input type="checkbox"/>	N/A	
		Short with casters	R895	\$270 <input type="checkbox"/>	N/A	
		Long with casters	R896	\$295 <input type="checkbox"/>	N/A	
* Pads	CHOOSE	Tan	R875	\$235 <input type="checkbox"/>	R875	\$235 <input type="checkbox"/>
		Blue	R876	\$235 <input type="checkbox"/>	R876	\$235 <input type="checkbox"/>
		Red	R877	\$235 <input type="checkbox"/>	R877	\$235 <input type="checkbox"/>
		Pink	R977	\$235 <input type="checkbox"/>	R977	\$235 <input type="checkbox"/>
		Purple	R978	\$235 <input type="checkbox"/>	R978	\$235 <input type="checkbox"/>
		Green	R979	\$235 <input type="checkbox"/>	R979	\$235 <input type="checkbox"/>
Headrest	CHOOSE	Contoured	R897	\$220 <input type="checkbox"/>	R897	\$220 <input type="checkbox"/>
		Flat	R898	\$205 <input type="checkbox"/>	R898	\$205 <input type="checkbox"/>
		Adjustable winged	R881	\$270 <input type="checkbox"/>	R881	\$270 <input type="checkbox"/>
Laterals (pair)	CHOOSE	Small	R828	\$180 <input type="checkbox"/>	R828	\$180 <input type="checkbox"/>
		Large	R838	\$195 <input type="checkbox"/>	R838	\$195 <input type="checkbox"/>
		Large with strap	R839	\$235 <input type="checkbox"/>	R839	\$235 <input type="checkbox"/>
Push handles (recommended for all chairs with wheels)		R878	\$140 <input type="checkbox"/>	Included		
Tray (requires armrests)		R669	\$260 <input type="checkbox"/>	R669	\$260 <input type="checkbox"/>	
Handhold		R649	\$22 <input type="checkbox"/>	R649	\$22 <input type="checkbox"/>	
Chest strap (wide)		R807	\$85 <input type="checkbox"/>	R807	\$85 <input type="checkbox"/>	
Butterfly harness	CHOOSE	Regular	R889	\$135 <input type="checkbox"/>	R889	\$135 <input type="checkbox"/>
		Slim cut	R826	\$135 <input type="checkbox"/>	R826	\$135 <input type="checkbox"/>
Thigh belt		R813	\$85 <input type="checkbox"/>	R813	\$85 <input type="checkbox"/>	
Pelvic harness		R846	\$160 <input type="checkbox"/>	R846	\$160 <input type="checkbox"/>	
Hip guides (pair)		R867	\$190 <input type="checkbox"/>	R867	\$190 <input type="checkbox"/>	
Abductor	CHOOSE	Small	R884	\$115 <input type="checkbox"/>	R884	\$115 <input type="checkbox"/>
		Large	R891	\$135 <input type="checkbox"/>	R891	\$135 <input type="checkbox"/>
Adductors (pair)		R892	\$180 <input type="checkbox"/>	R892	\$180 <input type="checkbox"/>	
Leg prompt		R859	\$215 <input type="checkbox"/>	R859	\$215 <input type="checkbox"/>	
Footboard (recommended for chairs with long legs)		R857	\$445 <input type="checkbox"/>	Included		
Footboard lift		R818	\$220 <input type="checkbox"/>	R818	\$220 <input type="checkbox"/>	
Ankle straps (pair)		R855	\$86 <input type="checkbox"/>	R855	\$86 <input type="checkbox"/>	
Sandals (pair)		R964	\$98 <input type="checkbox"/>	R964	\$98 <input type="checkbox"/>	
Wedges (pair, requires sandals)		R699	\$67 <input type="checkbox"/>	R699	\$67 <input type="checkbox"/>	
Lumbar and seat support kit		R888	\$80 <input type="checkbox"/>	R888	\$80 <input type="checkbox"/>	
Backrest filler pad		R856	\$115 <input type="checkbox"/>	R856	\$115 <input type="checkbox"/>	
Whitmyer headrest adapter		R804	\$50 <input type="checkbox"/>	R804	\$50 <input type="checkbox"/>	
Total		R860	\$ _____	R870	\$ _____	

* Required

Room to grow

Large chair at its smallest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 9½"
- 3 Seat depth 15"
- 4 Seat width 11"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 15" (with footboard lift 11")



Large chair at its largest dimensions


- 1 Backrest height 24"
- 2 Distance between laterals 14"
- 3 Seat depth 20"
- 4 Seat width 17"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 20"




Large standard base
R860 Rifton Activity Chair



Large hi/lo base
R870 Rifton Activity Chair

User dimensions (inches)	
Height	50–74
 <p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	
	R860 large Standard base
Frame width	short legs: 26¾ long legs: 26¾
	short legs w/ casters: 27¾ long legs w/ casters: 27¾
Seat height above floor	short legs: 16½–20½ long legs: 20–24
	short legs w/ casters: 16½–20½ long legs w/ casters: 20–24
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Seat height above footboard with footboard lift	15–20 11–16
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lbs)	250

User dimensions (inches)	
Height	50–74
 <p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	
	R870 large Hi/lo base
Frame width	29¾
Seat height above floor	16–24
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Seat height above footboard with footboard lift	15–20 11–16
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lbs)	250