

# 2021 order form for large

Use dimension chart to select appropriate size.



		Standard R860 large		Hi/lo R870 large		
* <b>Seat and back</b> includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R861	\$660 <input type="checkbox"/>	R871	\$1,260 <input type="checkbox"/>
		Adjustable w/ back spring	R862	\$765 <input type="checkbox"/>	R872	\$1,365 <input type="checkbox"/>
* <b>Base</b>	CHOOSE	Without spring	R873	\$405 <input type="checkbox"/>	R883	\$2,730 <input type="checkbox"/>
		With spring	R874	\$485 <input type="checkbox"/>	N/A	
* <b>Arm supports</b> (pair)	CHOOSE	Armrests	R868	\$235 <input type="checkbox"/>	R868	\$235 <input type="checkbox"/>
		Forearm prompts	R879	\$610 <input type="checkbox"/>	R879	\$610 <input type="checkbox"/>
* <b>Adjustable legs</b>	CHOOSE	Short	R808	\$65 <input type="checkbox"/>	N/A	
		Long	R894	\$75 <input type="checkbox"/>	N/A	
		Short with casters	R895	\$285 <input type="checkbox"/>	N/A	
		Long with casters	R896	\$310 <input type="checkbox"/>	N/A	
* <b>Pads</b>	CHOOSE	Tan	R875	\$245 <input type="checkbox"/>	R875	\$245 <input type="checkbox"/>
		Blue	R876	\$245 <input type="checkbox"/>	R876	\$245 <input type="checkbox"/>
		Red	R877	\$245 <input type="checkbox"/>	R877	\$245 <input type="checkbox"/>
		Pink	R977	\$245 <input type="checkbox"/>	R977	\$245 <input type="checkbox"/>
		Purple	R978	\$245 <input type="checkbox"/>	R978	\$245 <input type="checkbox"/>
		Green	R979	\$245 <input type="checkbox"/>	R979	\$245 <input type="checkbox"/>
<b>Headrest</b>	CHOOSE	Contoured	R897	\$230 <input type="checkbox"/>	R897	\$230 <input type="checkbox"/>
		Flat	R898	\$215 <input type="checkbox"/>	R898	\$215 <input type="checkbox"/>
		Adjustable winged	R881	\$285 <input type="checkbox"/>	R881	\$285 <input type="checkbox"/>
<b>Laterals</b> (pair)	CHOOSE	Small	R828	\$190 <input type="checkbox"/>	R828	\$190 <input type="checkbox"/>
		Large	R838	\$205 <input type="checkbox"/>	R838	\$205 <input type="checkbox"/>
		Large with strap	R839	\$245 <input type="checkbox"/>	R839	\$245 <input type="checkbox"/>
<b>Push handles</b> (recommended for all chairs with wheels)			R878	\$145 <input type="checkbox"/>	Included	
<b>Tray</b> (requires armrests)			R669	\$275 <input type="checkbox"/>	R669	\$275 <input type="checkbox"/>
<b>Handhold</b>			R649	\$25 <input type="checkbox"/>	R649	\$25 <input type="checkbox"/>
<b>Chest strap</b> (wide)			R807	\$90 <input type="checkbox"/>	R807	\$90 <input type="checkbox"/>
<b>Butterfly harness</b>	CHOOSE	Regular	R889	\$140 <input type="checkbox"/>	R889	\$140 <input type="checkbox"/>
		Slim cut	R826	\$140 <input type="checkbox"/>	R826	\$140 <input type="checkbox"/>
<b>Thigh belt</b>			R813	\$90 <input type="checkbox"/>	R813	\$90 <input type="checkbox"/>
<b>Pelvic harness</b>			R846	\$170 <input type="checkbox"/>	R846	\$170 <input type="checkbox"/>
<b>Hip guides</b> (pair)			R867	\$200 <input type="checkbox"/>	R867	\$200 <input type="checkbox"/>
<b>Abductor</b>	CHOOSE	Small	R884	\$120 <input type="checkbox"/>	R884	\$120 <input type="checkbox"/>
		Large	R891	\$140 <input type="checkbox"/>	R891	\$140 <input type="checkbox"/>
<b>Adductors</b> (pair)			R892	\$190 <input type="checkbox"/>	R892	\$190 <input type="checkbox"/>
<b>Leg prompt</b>			R859	\$225 <input type="checkbox"/>	R859	\$225 <input type="checkbox"/>
<b>Footboard</b> (recommended for chairs with long legs)			R857	\$470 <input type="checkbox"/>	Included	
<b>Ankle straps</b> (pair)			R855	\$90 <input type="checkbox"/>	R855	\$90 <input type="checkbox"/>
<b>Sandals</b> (pair)			R964	\$105 <input type="checkbox"/>	R964	\$105 <input type="checkbox"/>
<b>Wedges</b> (pair, require sandals)			R699	\$70 <input type="checkbox"/>	R699	\$70 <input type="checkbox"/>
<b>Lumbar and seat support kit</b>			R888	\$85 <input type="checkbox"/>	R888	\$85 <input type="checkbox"/>
<b>Backrest filler pad</b>			R856	\$120 <input type="checkbox"/>	R856	\$120 <input type="checkbox"/>
<b>Whitmyer headrest adapter</b>			R804	\$55 <input type="checkbox"/>	R804	\$55 <input type="checkbox"/>
<b>Total</b>			<b>R860</b>	<b>\$ _____</b>	<b>R870</b>	<b>\$ _____</b>

\* Required

# Room to grow

## Large chair at its smallest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 9½"
- 3 Seat depth 15"
- 4 Seat width 11"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 11¼"




## Large chair at its largest dimensions

- 1 Backrest height 24"
- 2 Distance between laterals 14"
- 3 Seat depth 20"
- 4 Seat width 17"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 21¼"




**Large standard base**  
R860 Rifton Activity Chair

User dimensions (inches)	
Height	50–74
 <p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	
	R860 large Standard base
Frame width	short legs: 26¾ long legs: 26¾
	short legs w/ casters: 27¾ long legs w/ casters: 27¾
Seat height above floor	short legs: 16½–20½ long legs: 20–24
	short legs w/ casters: 16½–20½ long legs w/ casters: 20–24
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lb)	250



**Large hi/lo base**  
R870 Rifton Activity Chair

User dimensions (inches)	
Height	50–74
 <p><b>Key user dimension: height</b> The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	
	R870 large Hi/lo base
Frame width	29¾
Seat height above floor	16–24
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	11¼–21¼
Seat width with hip guides without hip guides	11–14 17
Seat depth	15–20
Armrest height above seat	7–10½
Distance between laterals	9½–14
Backrest height	19–24
Headrest height above seat	19½–29½
Max. working load (lb)	250