

Medium Activity Chair 2018



			Standard R840 medium		Hi/lo R850 medium		
Required components	Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE Adjustable	R841	\$530 <input type="checkbox"/>	R851	\$1,035 <input type="checkbox"/>	
		Adjustable w/ back spring	R842	\$605 <input type="checkbox"/>	R852	\$1,135 <input type="checkbox"/>	
	Base	CHOOSE Without spring	R843	\$310 <input type="checkbox"/>	R853	\$2,225 <input type="checkbox"/>	
		With spring	R844	\$370 <input type="checkbox"/>	N/A		
	Arm supports (pair)	CHOOSE Armrests	R868	\$215 <input type="checkbox"/>	R868	\$215 <input type="checkbox"/>	
		Forearm prompts	R869	\$540 <input type="checkbox"/>	R869	\$540 <input type="checkbox"/>	
	Adjustable legs	CHOOSE	Short	R893	\$48 <input type="checkbox"/>	N/A	
			Long	R894	\$70 <input type="checkbox"/>	N/A	
			Short with casters	R895	\$260 <input type="checkbox"/>	N/A	
			Long with casters	R896	\$285 <input type="checkbox"/>	N/A	
	Pads	CHOOSE	Tan	R835	\$215 <input type="checkbox"/>	R835	\$215 <input type="checkbox"/>
			Blue	R836	\$215 <input type="checkbox"/>	R836	\$215 <input type="checkbox"/>
			Red	R837	\$215 <input type="checkbox"/>	R837	\$215 <input type="checkbox"/>
			Pink	R974	\$215 <input type="checkbox"/>	R974	\$215 <input type="checkbox"/>
Purple			R975	\$215 <input type="checkbox"/>	R975	\$215 <input type="checkbox"/>	
Green			R976	\$215 <input type="checkbox"/>	R976	\$215 <input type="checkbox"/>	
Optional accessories	Headrest	CHOOSE Contoured	R897	\$210 <input type="checkbox"/>	R897	\$210 <input type="checkbox"/>	
		Flat	R898	\$195 <input type="checkbox"/>	R898	\$195 <input type="checkbox"/>	
		Adjustable winged	R881	\$260 <input type="checkbox"/>	R881	\$260 <input type="checkbox"/>	
	Laterals (pair)	CHOOSE	Small	R828	\$175 <input type="checkbox"/>	R828	\$175 <input type="checkbox"/>
			Small with strap	R829	\$215 <input type="checkbox"/>	R829	\$215 <input type="checkbox"/>
			Large	R838	\$190 <input type="checkbox"/>	R838	\$190 <input type="checkbox"/>
	Push handles (recommended for all chairs with wheels)			R878	\$135 <input type="checkbox"/>	Included	
	Tray (requires armrests)			R668	\$240 <input type="checkbox"/>	R668	\$240 <input type="checkbox"/>
	Handhold			R649	\$21 <input type="checkbox"/>	R649	\$21 <input type="checkbox"/>
	Chest strap (wide)			R863	\$70 <input type="checkbox"/>	R863	\$70 <input type="checkbox"/>
	Butterfly harness	CHOOSE	Regular	R864	\$110 <input type="checkbox"/>	R864	\$110 <input type="checkbox"/>
			Slim cut	R882	\$110 <input type="checkbox"/>	R882	\$110 <input type="checkbox"/>
	Thigh belt			R865	\$70 <input type="checkbox"/>	R865	\$70 <input type="checkbox"/>
	Pelvic harness			R866	\$140 <input type="checkbox"/>	R866	\$140 <input type="checkbox"/>
	Hip guides (pair)			R867	\$185 <input type="checkbox"/>	R867	\$185 <input type="checkbox"/>
	Abductor	CHOOSE	Small	R884	\$110 <input type="checkbox"/>	R884	\$110 <input type="checkbox"/>
			Large	R891	\$130 <input type="checkbox"/>	R891	\$130 <input type="checkbox"/>
	Adductors (pair)			R892	\$175 <input type="checkbox"/>	R892	\$175 <input type="checkbox"/>
	Leg prompt			R858	\$195 <input type="checkbox"/>	R858	\$195 <input type="checkbox"/>
	Footboard (recommended for chairs with long legs)			R848	\$380 <input type="checkbox"/>	Included	
	Footboard lift			R806	\$205 <input type="checkbox"/>	R806	\$205 <input type="checkbox"/>
	Ankle straps (pair)			R899	\$70 <input type="checkbox"/>	R899	\$70 <input type="checkbox"/>
	Sandals (pair)			R954	\$83 <input type="checkbox"/>	R954	\$83 <input type="checkbox"/>
	Wedges (pair, requires sandals)			R698	\$59 <input type="checkbox"/>	R698	\$59 <input type="checkbox"/>
	Lumbar and seat support kit			R887	\$70 <input type="checkbox"/>	R887	\$70 <input type="checkbox"/>
	Backrest filler pad			R849	\$110 <input type="checkbox"/>	R849	\$110 <input type="checkbox"/>
	Whitmyer headrest adapter			R804	\$40 <input type="checkbox"/>	R804	\$40 <input type="checkbox"/>
Total			R840	\$ _____	R850	\$ _____	

Room to grow

Medium chair at its smallest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 6½"
- 3 Seat depth 11"
- 4 Seat width 8½"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 12"
(with footboard lift 8")



Medium chair at its largest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 12"
- 3 Seat depth 16"
- 4 Seat width 14"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 16½"



Medium standard base
R840 Rifton Activity Chair

User dimensions (inches)	
Height	40–56
<p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	
Frame width	short legs: 23¼
	long legs: 24½
	short legs w/ casters: 24½
	long legs w/ casters: 25½
Seat height above floor	short legs: 12½–16½
	long legs: 19–23
	short legs w/ casters: 16–20
	long legs w/ casters: 19–23
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Seat height above footboard with footboard lift	12–16½
	8–12½
Seat width with hip guides without hip guides	8½–11½
	14
Seat depth	11–16
Armrest height above seat	7–10½
Distance between laterals	6½–12
Backrest height	15½–19
Headrest height above seat	17½–24½
Max. working load (lbs)	150



Medium hi/lo base
R850 Rifton Activity Chair

User dimensions (inches)	
Height	40–56
<p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	
Frame width	27½
Seat height above floor	12–23
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Seat height above footboard with footboard lift	12–16½
	8–12½
Seat width with hip guides without hip guides	8½–11½
	14
Seat depth	11–16
Armrest height above seat	7–10½
Distance between laterals	6½–12
Backrest height	15½–19
Headrest height above seat	17½–24½
Max. working load (lbs)	150