

2019 order form for medium

Use dimension chart to select appropriate size.



		Standard R840 medium		Hi/lo R850 medium		
* Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R841	\$550 <input type="checkbox"/>	R851	\$1,150 <input type="checkbox"/>
		Adjustable w/ back spring	R842	\$630 <input type="checkbox"/>	R852	\$1,250 <input type="checkbox"/>
* Base	CHOOSE	Without spring	R843	\$320 <input type="checkbox"/>	R853	\$2,315 <input type="checkbox"/>
		With spring	R844	\$385 <input type="checkbox"/>	N/A	
* Arm supports (pair)	CHOOSE	Armrests	R868	\$225 <input type="checkbox"/>	R868	\$225 <input type="checkbox"/>
		Forearm prompts	R869	\$560 <input type="checkbox"/>	R869	\$560 <input type="checkbox"/>
* Adjustable legs	CHOOSE	Short	R893	\$50 <input type="checkbox"/>	N/A	
		Long	R894	\$73 <input type="checkbox"/>	N/A	
		Short with casters	R895	\$270 <input type="checkbox"/>	N/A	
		Long with casters	R896	\$295 <input type="checkbox"/>	N/A	
* Pads	CHOOSE	Tan	R835	\$225 <input type="checkbox"/>	R835	\$225 <input type="checkbox"/>
		Blue	R836	\$225 <input type="checkbox"/>	R836	\$225 <input type="checkbox"/>
		Red	R837	\$225 <input type="checkbox"/>	R837	\$225 <input type="checkbox"/>
		Pink	R974	\$225 <input type="checkbox"/>	R974	\$225 <input type="checkbox"/>
		Purple	R975	\$225 <input type="checkbox"/>	R975	\$225 <input type="checkbox"/>
		Green	R976	\$225 <input type="checkbox"/>	R976	\$225 <input type="checkbox"/>
Headrest	CHOOSE	Contoured	R897	\$220 <input type="checkbox"/>	R897	\$220 <input type="checkbox"/>
		Flat	R898	\$205 <input type="checkbox"/>	R898	\$205 <input type="checkbox"/>
		Adjustable winged	R881	\$270 <input type="checkbox"/>	R881	\$270 <input type="checkbox"/>
Laterals (pair)	CHOOSE	Small	R828	\$180 <input type="checkbox"/>	R828	\$180 <input type="checkbox"/>
		Small with strap	R829	\$225 <input type="checkbox"/>	R829	\$225 <input type="checkbox"/>
		Large	R838	\$195 <input type="checkbox"/>	R838	\$195 <input type="checkbox"/>
Push handles (recommended for all chairs with wheels)			R878	\$140 <input type="checkbox"/>	Included	
Tray (requires armrests)			R668	\$250 <input type="checkbox"/>	R668	\$250 <input type="checkbox"/>
Handhold			R649	\$22 <input type="checkbox"/>	R649	\$22 <input type="checkbox"/>
Chest strap (wide)			R863	\$80 <input type="checkbox"/>	R863	\$80 <input type="checkbox"/>
Butterfly harness	CHOOSE	Regular	R864	\$115 <input type="checkbox"/>	R864	\$115 <input type="checkbox"/>
		Slim cut	R882	\$115 <input type="checkbox"/>	R882	\$115 <input type="checkbox"/>
Thigh belt			R865	\$80 <input type="checkbox"/>	R865	\$80 <input type="checkbox"/>
Pelvic harness			R866	\$145 <input type="checkbox"/>	R866	\$145 <input type="checkbox"/>
Hip guides (pair)			R867	\$190 <input type="checkbox"/>	R867	\$190 <input type="checkbox"/>
Abductor	CHOOSE	Small	R884	\$115 <input type="checkbox"/>	R884	\$115 <input type="checkbox"/>
		Large	R891	\$135 <input type="checkbox"/>	R891	\$135 <input type="checkbox"/>
Adductors (pair)			R892	\$180 <input type="checkbox"/>	R892	\$180 <input type="checkbox"/>
Leg prompt			R858	\$205 <input type="checkbox"/>	R858	\$205 <input type="checkbox"/>
Footboard (recommended for chairs with long legs)			R848	\$435 <input type="checkbox"/>	Included	
Footboard lift			R806	\$215 <input type="checkbox"/>	R806	\$215 <input type="checkbox"/>
Ankle straps (pair)			R899	\$80 <input type="checkbox"/>	R899	\$80 <input type="checkbox"/>
Sandals (pair)			R954	\$86 <input type="checkbox"/>	R954	\$86 <input type="checkbox"/>
Wedges (pair, requires sandals)			R698	\$61 <input type="checkbox"/>	R698	\$61 <input type="checkbox"/>
Lumbar and seat support kit			R887	\$73 <input type="checkbox"/>	R887	\$73 <input type="checkbox"/>
Backrest filler pad			R849	\$115 <input type="checkbox"/>	R849	\$115 <input type="checkbox"/>
Whitmyer headrest adapter			R804	\$50 <input type="checkbox"/>	R804	\$50 <input type="checkbox"/>
Total			R840	\$ _____	R850	\$ _____

* Required

Room to grow

Medium chair at its smallest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 6½"
- 3 Seat depth 11"
- 4 Seat width 8½"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 12"
(with footboard lift 8")



Medium chair at its largest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 12"
- 3 Seat depth 16"
- 4 Seat width 14"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 16½"



Medium standard base
R840 Rifton Activity Chair

User dimensions (inches)	
Height	40–56
	Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.
	Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.
Item dimensions (inches)	R840 medium Standard base
Frame width	short legs: 23¼
	long legs: 24½
	short legs w/ casters: 24½
	long legs w/ casters: 25½
Seat height above floor	short legs: 12½–16½
	long legs: 19–23
	short legs w/ casters: 16–20
	long legs w/ casters: 19–23
Seat angle (tilt-in-space)–degrees	15° forward, 15° back
Backrest angle–degrees	5° forward, 20° back
Seat height above footboard with footboard lift	12–16½
	8–12½
Seat width with hip guides without hip guides	8½–11½
	14
Seat depth	11–16
Armrest height above seat	7–10½
Distance between laterals	6½–12
Backrest height	15½–19
Headrest height above seat	17½–24½
Max. working load (lbs)	150



Medium hi/lo base
R850 Rifton Activity Chair

User dimensions (inches)	
Height	40–56
	Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.
	Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.
Item dimensions (inches)	R850 medium Hi/lo base
Frame width	27½
Seat height above floor	12–23
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Seat height above footboard with footboard lift	12–16½
	8–12½
Seat width with hip guides without hip guides	8½–11½
	14
Seat depth	11–16
Armrest height above seat	7–10½
Distance between laterals	6½–12
Backrest height	15½–19
Headrest height above seat	17½–24½
Max. working load (lbs)	150