

2020 order form for medium

Use dimension chart to select appropriate size.



		Standard R840 medium		Hi/lo R850 medium		
* Seat and back includes seatbelt (Hi/lo includes footboard and push handles)	CHOOSE	Adjustable	R841	\$580 <input type="checkbox"/>	R851	\$1,210 <input type="checkbox"/>
		Adjustable w/ back spring	R842	\$660 <input type="checkbox"/>	R852	\$1,315 <input type="checkbox"/>
* Base	CHOOSE	Without spring	R843	\$335 <input type="checkbox"/>	R853	\$2,430 <input type="checkbox"/>
		With spring	R844	\$405 <input type="checkbox"/>	N/A	
* Arm supports (pair)	CHOOSE	Armrests	R868	\$235 <input type="checkbox"/>	R868	\$235 <input type="checkbox"/>
		Forearm prompts	R869	\$590 <input type="checkbox"/>	R869	\$590 <input type="checkbox"/>
* Adjustable legs	CHOOSE	Short	R893	\$55 <input type="checkbox"/>	N/A	
		Long	R894	\$75 <input type="checkbox"/>	N/A	
		Short with casters	R895	\$285 <input type="checkbox"/>	N/A	
		Long with casters	R896	\$310 <input type="checkbox"/>	N/A	
* Pads	CHOOSE	Tan	R835	\$235 <input type="checkbox"/>	R835	\$235 <input type="checkbox"/>
		Blue	R836	\$235 <input type="checkbox"/>	R836	\$235 <input type="checkbox"/>
		Red	R837	\$235 <input type="checkbox"/>	R837	\$235 <input type="checkbox"/>
		Pink	R974	\$235 <input type="checkbox"/>	R974	\$235 <input type="checkbox"/>
		Purple	R975	\$235 <input type="checkbox"/>	R975	\$235 <input type="checkbox"/>
		Green	R976	\$235 <input type="checkbox"/>	R976	\$235 <input type="checkbox"/>
Headrest	CHOOSE	Contoured	R897	\$230 <input type="checkbox"/>	R897	\$230 <input type="checkbox"/>
		Flat	R898	\$215 <input type="checkbox"/>	R898	\$215 <input type="checkbox"/>
		Adjustable winged	R881	\$285 <input type="checkbox"/>	R881	\$285 <input type="checkbox"/>
Laterals (pair)	CHOOSE	Small	R828	\$190 <input type="checkbox"/>	R828	\$190 <input type="checkbox"/>
		Small with strap	R829	\$235 <input type="checkbox"/>	R829	\$235 <input type="checkbox"/>
		Large	R838	\$205 <input type="checkbox"/>	R838	\$205 <input type="checkbox"/>
Push handles (recommended for all chairs with wheels)		R878	\$145 <input type="checkbox"/>	Included		
Tray (requires armrests)		R668	\$265 <input type="checkbox"/>	R668	\$265 <input type="checkbox"/>	
Handhold		R649	\$25 <input type="checkbox"/>	R649	\$25 <input type="checkbox"/>	
Chest strap (wide)		R863	\$85 <input type="checkbox"/>	R863	\$85 <input type="checkbox"/>	
Butterfly harness	CHOOSE	Regular	R864	\$120 <input type="checkbox"/>	R864	\$120 <input type="checkbox"/>
		Slim cut	R882	\$120 <input type="checkbox"/>	R882	\$120 <input type="checkbox"/>
Thigh belt		R865	\$85 <input type="checkbox"/>	R865	\$85 <input type="checkbox"/>	
Pelvic harness		R866	\$150 <input type="checkbox"/>	R866	\$150 <input type="checkbox"/>	
Hip guides (pair)		R867	\$200 <input type="checkbox"/>	R867	\$200 <input type="checkbox"/>	
Abductor	CHOOSE	Small	R884	\$120 <input type="checkbox"/>	R884	\$120 <input type="checkbox"/>
		Large	R891	\$140 <input type="checkbox"/>	R891	\$140 <input type="checkbox"/>
Adductors (pair)		R892	\$190 <input type="checkbox"/>	R892	\$190 <input type="checkbox"/>	
Leg prompt		R858	\$215 <input type="checkbox"/>	R858	\$215 <input type="checkbox"/>	
Footboard (recommended for chairs with long legs)		R848	\$455 <input type="checkbox"/>	Included		
Ankle straps (pair)		R899	\$85 <input type="checkbox"/>	R899	\$85 <input type="checkbox"/>	
Sandals (pair)		R954	\$90 <input type="checkbox"/>	R954	\$90 <input type="checkbox"/>	
Wedges (pair, requires sandals)		R698	\$65 <input type="checkbox"/>	R698	\$65 <input type="checkbox"/>	
Lumbar and seat support kit		R887	\$75 <input type="checkbox"/>	R887	\$75 <input type="checkbox"/>	
Backrest filler pad		R849	\$120 <input type="checkbox"/>	R849	\$120 <input type="checkbox"/>	
Whitmyer headrest adapter		R804	\$55 <input type="checkbox"/>	R804	\$55 <input type="checkbox"/>	
Total		R840	\$ _____	R850	\$ _____	

* Required

Room to grow

Medium chair at its smallest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 6½"
- 3 Seat depth 11"
- 4 Seat width 8½"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 7¾"




Medium chair at its largest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 12"
- 3 Seat depth 16"
- 4 Seat width 14"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 15¼"




Medium standard base
R840 Rifton Activity Chair

User dimensions (inches)	
Height	40–56
 <p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	
Frame width	short legs: 23¾
	long legs: 24½
Seat height above floor	short legs w/ casters: 24½
	long legs w/ casters: 25½
Seat height above footboard	short legs: 12½–16½
	long legs: 19–23
Seat angle (tilt-in-space)–degrees	short legs w/ casters: 16–20
	long legs w/ casters: 19–23
Backrest angle–degrees	15° forward, 15° back
Footboard knee angle–degrees	5° forward, 20° back
Seat height above footboard	45°–110°
Seat width with hip guides	7¾–15¼
without hip guides	8½–11½
Seat depth	14
Armrest height above seat	11–16
Distance between laterals	7–10½
Backrest height	6½–12
Headrest height above seat	15½–19
Max. working load (lbs)	17½–24½
	150



Medium hi/lo base
R850 Rifton Activity Chair

User dimensions (inches)	
Height	40–56
 <p>Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p>Important: Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p>	
Item dimensions (inches)	
Frame width	27½
Seat height above floor	12–23
Seat angle (tilt-in-space)–degrees	15° forward, 25° back
Backrest angle–degrees	5° forward, 20° back
Footboard knee angle–degrees	45°–110°
Seat height above footboard	7¾–15¼
Seat width with hip guides	8½–11½
without hip guides	14
Seat depth	11–16
Armrest height above seat	7–10½
Distance between laterals	6½–12
Backrest height	15½–19
Headrest height above seat	17½–24½
Max. working load (lbs)	150