

# 2021 order form for medium

Use dimension chart to select appropriate size.



|  |        | Standard<br>R840 medium   |             | Hi/lo<br>R850 medium           |             |                                  |
|--|--------|---------------------------|-------------|--------------------------------|-------------|----------------------------------|
| * <b>Seat and back</b> includes seatbelt<br>(Hi/lo includes footboard and<br>push handles) | CHOOSE | Adjustable                | R841        | \$580 <input type="checkbox"/> | R851        | \$1,210 <input type="checkbox"/> |
|  |        | Adjustable w/ back spring | R842        | \$660 <input type="checkbox"/> | R852        | \$1,315 <input type="checkbox"/> |
| * <b>Base</b>  | CHOOSE | Without spring            | R843        | \$335 <input type="checkbox"/> | R853        | \$2,430 <input type="checkbox"/> |
|  |        | With spring               | R844        | \$405 <input type="checkbox"/> | N/A         |                                  |
| * <b>Arm supports</b> (pair)   | CHOOSE | Armrests                  | R868        | \$235 <input type="checkbox"/> | R868        | \$235 <input type="checkbox"/>   |
|  |        | Forearm prompts           | R869        | \$590 <input type="checkbox"/> | R869        | \$590 <input type="checkbox"/>   |
| * <b>Adjustable legs</b>   | CHOOSE | Short                     | R893        | \$55 <input type="checkbox"/>  | N/A         |                                  |
|  |        | Long                      | R894        | \$75 <input type="checkbox"/>  | N/A         |                                  |
|  |        | Short with casters        | R895        | \$285 <input type="checkbox"/> | N/A         |                                  |
|  |        | Long with casters         | R896        | \$310 <input type="checkbox"/> | N/A         |                                  |
| * <b>Pads</b>  | CHOOSE | Tan                       | R835        | \$235 <input type="checkbox"/> | R835        | \$235 <input type="checkbox"/>   |
|  |        | Blue                      | R836        | \$235 <input type="checkbox"/> | R836        | \$235 <input type="checkbox"/>   |
|  |        | Red                       | R837        | \$235 <input type="checkbox"/> | R837        | \$235 <input type="checkbox"/>   |
|  |        | Pink                      | R974        | \$235 <input type="checkbox"/> | R974        | \$235 <input type="checkbox"/>   |
|  |        | Purple                    | R975        | \$235 <input type="checkbox"/> | R975        | \$235 <input type="checkbox"/>   |
|  |        | Green                     | R976        | \$235 <input type="checkbox"/> | R976        | \$235 <input type="checkbox"/>   |
| <b>Headrest</b>  | CHOOSE | Contoured                 | R897        | \$230 <input type="checkbox"/> | R897        | \$230 <input type="checkbox"/>   |
|  |        | Flat                      | R898        | \$215 <input type="checkbox"/> | R898        | \$215 <input type="checkbox"/>   |
|  |        | Adjustable winged         | R881        | \$285 <input type="checkbox"/> | R881        | \$285 <input type="checkbox"/>   |
| <b>Laterals</b> (pair)   | CHOOSE | Small                     | R828        | \$190 <input type="checkbox"/> | R828        | \$190 <input type="checkbox"/>   |
|  |        | Small with strap          | R829        | \$235 <input type="checkbox"/> | R829        | \$235 <input type="checkbox"/>   |
|  |        | Large                     | R838        | \$205 <input type="checkbox"/> | R838        | \$205 <input type="checkbox"/>   |
| <b>Push handles</b> (recommended for all chairs with wheels)                               |        |                           | R878        | \$145 <input type="checkbox"/> | Included    |                                  |
| <b>Tray</b> (requires armrests)  |        |                           | R668        | \$265 <input type="checkbox"/> | R668        | \$265 <input type="checkbox"/>   |
| <b>Handhold</b>  |        |                           | R649        | \$25 <input type="checkbox"/>  | R649        | \$25 <input type="checkbox"/>    |
| <b>Chest strap</b> (wide)  |        |                           | R863        | \$85 <input type="checkbox"/>  | R863        | \$85 <input type="checkbox"/>    |
| <b>Butterfly harness</b>   | CHOOSE | Regular                   | R864        | \$120 <input type="checkbox"/> | R864        | \$120 <input type="checkbox"/>   |
|  |        | Slim cut                  | R882        | \$120 <input type="checkbox"/> | R882        | \$120 <input type="checkbox"/>   |
| <b>Thigh belt</b>  |        |                           | R865        | \$85 <input type="checkbox"/>  | R865        | \$85 <input type="checkbox"/>    |
| <b>Pelvic harness</b>  |        |                           | R866        | \$150 <input type="checkbox"/> | R866        | \$150 <input type="checkbox"/>   |
| <b>Hip guides</b> (pair)   |        |                           | R867        | \$200 <input type="checkbox"/> | R867        | \$200 <input type="checkbox"/>   |
| <b>Abductor</b>  | CHOOSE | Small                     | R884        | \$120 <input type="checkbox"/> | R884        | \$120 <input type="checkbox"/>   |
|  |        | Large                     | R891        | \$140 <input type="checkbox"/> | R891        | \$140 <input type="checkbox"/>   |
| <b>Adductors</b> (pair)  |        |                           | R892        | \$190 <input type="checkbox"/> | R892        | \$190 <input type="checkbox"/>   |
| <b>Leg prompt</b>  |        |                           | R858        | \$215 <input type="checkbox"/> | R858        | \$215 <input type="checkbox"/>   |
| <b>Footboard</b> (recommended for chairs with long legs)                                   |        |                           | R848        | \$455 <input type="checkbox"/> | Included    |                                  |
| <b>Ankle straps</b> (pair)   |        |                           | R899        | \$85 <input type="checkbox"/>  | R899        | \$85 <input type="checkbox"/>    |
| <b>Sandals</b> (pair)  |        |                           | R954        | \$90 <input type="checkbox"/>  | R954        | \$90 <input type="checkbox"/>    |
| <b>Wedges</b> (pair, require sandals)  |        |                           | R698        | \$65 <input type="checkbox"/>  | R698        | \$65 <input type="checkbox"/>    |
| <b>Lumbar and seat support kit</b>   |        |                           | R887        | \$75 <input type="checkbox"/>  | R887        | \$75 <input type="checkbox"/>    |
| <b>Backrest filler pad</b>   |        |                           | R849        | \$120 <input type="checkbox"/> | R849        | \$120 <input type="checkbox"/>   |
| <b>Whitmyer headrest adapter</b>   |        |                           | R804        | \$55 <input type="checkbox"/>  | R804        | \$55 <input type="checkbox"/>    |
| <b>Total</b>   |        |                           | <b>R840</b> | <b>\$ _____</b>                | <b>R850</b> | <b>\$ _____</b>                  |

\* Required

# Room to grow

## Medium chair at its smallest dimensions

- 1 Backrest height 15½"
- 2 Distance between laterals 6½"
- 3 Seat depth 11"
- 4 Seat width 8½"
- 5 Armrest height above seat 7"
- 6 Seat height above footboard 7¾"




## Medium chair at its largest dimensions

- 1 Backrest height 19"
- 2 Distance between laterals 12"
- 3 Seat depth 16"
- 4 Seat width 14"
- 5 Armrest height above seat 10½"
- 6 Seat height above footboard 15¼"




**Medium standard base**  
R840 Rifton Activity Chair

| User dimensions (inches)  |                              |
|---|------------------------------|
| Height  | 40–56                        |
|  <p><b>Key user dimension: height</b><br/>The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p> |                              |
|   |                              |
| Item dimensions (inches)  |                              |
| Frame width   | short legs: 23¾              |
|   | long legs: 24½               |
| Seat height above floor   | short legs w/ casters: 24½   |
|   | long legs w/ casters: 25½    |
| Seat height above footboard   | short legs: 12½–16½          |
|   | long legs: 19–23             |
| Seat angle (tilt-in-space)–degrees  | short legs w/ casters: 16–20 |
|   | long legs w/ casters: 19–23  |
| Backrest angle–degrees  | 15° forward, 15° back        |
| Footboard knee angle–degrees  | 5° forward, 20° back         |
| Seat height above footboard   | 45°–110°                     |
| Seat width with hip guides  | 7¾–15¼                       |
| without hip guides  | 8½–11½                       |
| Seat depth  | 14                           |
| Armrest height above seat   | 11–16                        |
| Distance between laterals   | 7–10½                        |
| Backrest height   | 6½–12                        |
| Headrest height above seat  | 15½–19                       |
| Max. working load (lb)  | 17½–24½                      |
|   | 150                          |



**Medium hi/lo base**  
R850 Rifton Activity Chair

| User dimensions (inches)  |                       |
|---|-----------------------|
| Height  | 40–56                 |
|  <p><b>Key user dimension: height</b><br/>The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth.</p> <p><b>Important:</b> Make sure that seat width, depth and height are adequate for each individual user, and that the user's weight does not exceed the maximum weight recommended.</p> |                       |
|   |                       |
| Item dimensions (inches)  |                       |
| Frame width   | 27½                   |
| Seat height above floor   | 12–23                 |
| Seat angle (tilt-in-space)–degrees  | 15° forward, 25° back |
| Backrest angle–degrees  | 5° forward, 20° back  |
| Footboard knee angle–degrees  | 45°–110°              |
| Seat height above footboard   | 7¾–15¼                |
| Seat width with hip guides  | 8½–11½                |
| without hip guides  | 14                    |
| Seat depth  | 11–16                 |
| Armrest height above seat   | 7–10½                 |
| Distance between laterals   | 6½–12                 |
| Backrest height   | 15½–19                |
| Headrest height above seat  | 17½–24½               |
| Max. working load (lb)  | 150                   |