


# TRAM 2018

- 1. Choose** standard or low-base frame.
- 2. Select** a model with or without scale.
- 3. Choose** thigh straps. Then select additional accessories as your needs require.
- 4. Add** up the total cost.
- 5. Order** online in minutes at **Rifton.com** or see the back of this catalog for other ways to order.
- 6. Call** us any time at **800.571.8198** if you need assistance. We're always here to help.




		<b>K310 Standard</b> <small>(better for gait training)</small>		<b>K320 Low-base</b> <small>(fits under low beds)</small>				
<b>Required components</b>	<b>Frame</b> includes one battery, one battery charger, accessories tote and charger power cord (international power cords available on request)	CHOOSE	Frame with scale	K301	\$5,370 <input type="checkbox"/>	K321	\$5,370 <input type="checkbox"/>	
			Frame	K302	\$4,400 <input type="checkbox"/>	K322	\$4,400 <input type="checkbox"/>	
	<b>Thigh straps</b> (pair)	CHOOSE	Narrow	K308	\$68 <input type="checkbox"/>	K308	\$68 <input type="checkbox"/>	
			Wide	K338	\$68 <input type="checkbox"/>	K338	\$68 <input type="checkbox"/>	
<b>Optional accessories</b>	<b>New! Switch pole</b>			K666	\$205 <input type="checkbox"/>	K666	\$205 <input type="checkbox"/>	
	<b>Arm supports</b> (pair)	CHOOSE	Small arm prompts	K305	\$450 <input type="checkbox"/>	K305	\$450 <input type="checkbox"/>	
			Large arm prompts	K315	\$520 <input type="checkbox"/>	K315	\$520 <input type="checkbox"/>	
			Arm platforms	K397	\$335 <input type="checkbox"/>	K397	\$335 <input type="checkbox"/>	
			Arm platforms w/ handgrips	K396	\$385 <input type="checkbox"/>	K396	\$385 <input type="checkbox"/>	
	<b>Pelvic support</b>	CHOOSE	Small	K306	\$230 <input type="checkbox"/>	K306	\$230 <input type="checkbox"/>	
			Medium	K316	\$255 <input type="checkbox"/>	K316	\$255 <input type="checkbox"/>	
			Large	K326	\$270 <input type="checkbox"/>	K326	\$270 <input type="checkbox"/>	
	<b>New! Hip positioner</b>	CHOOSE	Small	K667	\$280 <input type="checkbox"/>	K667	\$280 <input type="checkbox"/>	
			Large	K668	\$310 <input type="checkbox"/>	K668	\$310 <input type="checkbox"/>	
	<b>New! Hip positioner pad</b>	CHOOSE	Small	K541	\$98 <input type="checkbox"/>	K541	\$98 <input type="checkbox"/>	
			Large	K542	\$130 <input type="checkbox"/>	K542	\$130 <input type="checkbox"/>	
	<b>Additional battery</b> (recommended for institutional use)				K307	\$125 <input type="checkbox"/>	K307	\$125 <input type="checkbox"/>
	<b>Additional thigh straps</b>	CHOOSE	Narrow single	K318	\$34 <input type="checkbox"/>	K318	\$34 <input type="checkbox"/>	
			Narrow pair	K308	\$68 <input type="checkbox"/>	K308	\$68 <input type="checkbox"/>	
Wide single			K328	\$34 <input type="checkbox"/>	K328	\$34 <input type="checkbox"/>		
Wide pair			K338	\$68 <input type="checkbox"/>	K338	\$68 <input type="checkbox"/>		
<b>Directional lock</b> (pair)				K319	\$170 <input type="checkbox"/>	K324	\$170 <input type="checkbox"/>	
<b>Removable back belt</b>				K327	\$53 <input type="checkbox"/>	K327	\$53 <input type="checkbox"/>	
<b>Total</b>				<b>K310</b>	<b>\$ _____</b>	<b>K320</b>	<b>\$ _____</b>	


# Dimensions of the TRAM

User dimensions (inches)		K310 and K320	
Girth		22 min – 60 max	
 <p><b>Key user dimension: girth</b></p> <p><b>Important:</b> User's weight or height must not exceed 350 lbs or 76 inches</p>			
Item dimensions (inches)	K310 TRAM	K320 TRAM	
Overall length	45½	Same	
Overall width	27½ – 40	Same	
Overall height	43½ – 58	Same	
TRAM weight (lbs)	70	Same	
Turning diameter	50	Same	
Floor to top of base leg for furniture clearance	6¾	4¾	
Min. user armpit height (when standing or ambulating)	30½	Same	
Max. user height (when standing or ambulating)	76	Same	
Max. working load (lbs)	350	Same	

## Pelvic support

User dimensions (inches) Indicated by the red arrow below			
<b>Key user dimension:</b> half of thigh circumference.			
K306 small	K316 med.	K326 large	
4 – 8	8 – 14	10 – 18	
			
Product dimensions (inches)			
K306 small	K316 med.	K326 large	
<b>a</b>	8½	11½	13½
<b>b</b>	15½	19¾	24
<b>c</b>	10	12½	15

## Hip positioner

Product dimensions (inches) – without pad		
		
	K513 small	K533 large
<b>a</b> Crotch length	2½	3½
<b>b</b> Width	1¼	1¾
<b>c</b> Length	9	11