

2019 order form


Use dimension chart to select appropriate size.




		K310 Standard (better for gait training)		K320 Low-base (fits under low beds)		
* Frame includes one battery, one battery charger, accessories tote and charger power cord (international power cords available on request)	CHOOSE	Frame with scale	K301	\$5,585 <input type="checkbox"/>	K321	\$5,585 <input type="checkbox"/>
		Frame	K302	\$4,575 <input type="checkbox"/>	K322	\$4,575 <input type="checkbox"/>
* Thigh straps (pair)	CHOOSE	Narrow	K308	\$71 <input type="checkbox"/>	K308	\$71 <input type="checkbox"/>
		Wide	K338	\$71 <input type="checkbox"/>	K338	\$71 <input type="checkbox"/>
Switch pole			K666	\$215 <input type="checkbox"/>	K666	\$215 <input type="checkbox"/>
Arm supports (pair)	CHOOSE	Small arm prompts	K305	\$450 <input type="checkbox"/>	K305	\$450 <input type="checkbox"/>
		Large arm prompts	K315	\$520 <input type="checkbox"/>	K315	\$520 <input type="checkbox"/>
		Arm platforms	K397	\$335 <input type="checkbox"/>	K397	\$335 <input type="checkbox"/>
		New! Arm platforms w/ handgrips	K396	\$385 <input type="checkbox"/>	K396	\$385 <input type="checkbox"/>
Pelvic support	CHOOSE	Small	K306	\$230 <input type="checkbox"/>	K306	\$230 <input type="checkbox"/>
		Medium	K316	\$255 <input type="checkbox"/>	K316	\$255 <input type="checkbox"/>
		Large	K326	\$270 <input type="checkbox"/>	K326	\$270 <input type="checkbox"/>
Hip positioner	CHOOSE	Small	K667	\$280 <input type="checkbox"/>	K667	\$280 <input type="checkbox"/>
		Large	K668	\$310 <input type="checkbox"/>	K668	\$310 <input type="checkbox"/>
Hip positioner pad	CHOOSE	Small	K541	\$98 <input type="checkbox"/>	K541	\$98 <input type="checkbox"/>
		Large	K542	\$130 <input type="checkbox"/>	K542	\$130 <input type="checkbox"/>
Additional battery (recommended for institutional use)			K307	\$130 <input type="checkbox"/>	K307	\$130 <input type="checkbox"/>
Additional thigh straps	CHOOSE	Narrow single	K318	\$35 <input type="checkbox"/>	K318	\$35 <input type="checkbox"/>
		Narrow pair	K308	\$71 <input type="checkbox"/>	K308	\$71 <input type="checkbox"/>
		Wide single	K328	\$35 <input type="checkbox"/>	K328	\$35 <input type="checkbox"/>
		Wide pair	K338	\$71 <input type="checkbox"/>	K338	\$71 <input type="checkbox"/>
Swivel lock (pair)			K319	\$177 <input type="checkbox"/>	K324	\$177 <input type="checkbox"/>
Removable back belt			K327	\$55 <input type="checkbox"/>	K327	\$55 <input type="checkbox"/>
Total			K310	\$ _____	K320	\$ _____

* Required

TRAM dimensions

User dimensions (inches)		K310 and K320	
Girth		22 min – 60 max	
 <p>Key user dimension: girth</p> <p>Important: User's weight or height must not exceed 350 lbs or 76 inches</p>			
Item dimensions (inches)	K310 TRAM	K320 TRAM	
Overall length	45½	Same	
Overall width	27½–40	Same	
Overall height	43½–58	Same	
TRAM weight (lbs)	70	Same	
Turning diameter	50	Same	
Floor to top of base leg for furniture clearance	6¾	4¼	
Min. user armpit height (when standing or ambulating)	30½	Same	
Max. user height (when standing or ambulating)	76	Same	
Max. working load (lbs)	350	Same	

Pelvic support

User dimensions (inches)			
Indicated by the red arrow below			
Key user dimension: half of thigh circumference measured at the groin.			
K306 small	K316 med.	K326 large	
6–10	8–14	10–18	
			
Product dimensions (inches)			
K306 small	K316 med.	K326 large	
a	8½	11½	13½
b	15½	19¾	24
c	10	12½	15

Hip positioner

Product dimensions (inches) – without pad		
	K667 small	K668 large
a Crotch length	2½	3½
b Width	1¼	1¾
c Length	9	11

