

# 2021 order form


Use dimension chart to select appropriate size.




		K310 Standard (better for gait training)		K320 Low-base (fits under low beds)	
* <b>Frame</b> includes one battery, one battery charger, accessories tote and charger power cord. (International power cords available on request.)	CHOOSE	Frame with scale	K301 \$5,860 <input type="checkbox"/>	K321 \$5,860 <input type="checkbox"/>	
		Frame	K302 \$4,800 <input type="checkbox"/>	K322 \$4,800 <input type="checkbox"/>	
* <b>Thigh straps</b> (pair)	CHOOSE	Narrow	K308 \$100 <input type="checkbox"/>	K308 \$100 <input type="checkbox"/>	
		Wide	K338 \$100 <input type="checkbox"/>	K338 \$100 <input type="checkbox"/>	
<b>Switch pole</b>		K666 \$225 <input type="checkbox"/>	K666 \$225 <input type="checkbox"/>		
<b>Arm supports</b> (pair)	CHOOSE	Small arm prompts	K305 \$475 <input type="checkbox"/>	K305 \$475 <input type="checkbox"/>	
		Large arm prompts	K315 \$545 <input type="checkbox"/>	K315 \$545 <input type="checkbox"/>	
		Arm platforms	K397 \$350 <input type="checkbox"/>	K397 \$350 <input type="checkbox"/>	
		Arm platforms w/ handgrips	K396 \$400 <input type="checkbox"/>	K396 \$400 <input type="checkbox"/>	
<b>Pelvic support</b>	CHOOSE	Small	K306 \$240 <input type="checkbox"/>	K306 \$240 <input type="checkbox"/>	
		Medium	K316 \$270 <input type="checkbox"/>	K316 \$270 <input type="checkbox"/>	
		Large	K326 \$285 <input type="checkbox"/>	K326 \$285 <input type="checkbox"/>	
<b>Hip positioner</b>	CHOOSE	Small	K667 \$295 <input type="checkbox"/>	K667 \$295 <input type="checkbox"/>	
		Large	K668 \$325 <input type="checkbox"/>	K668 \$325 <input type="checkbox"/>	
<b>Hip positioner pad</b>	CHOOSE	Small	K541 \$100 <input type="checkbox"/>	K541 \$100 <input type="checkbox"/>	
		Large	K542 \$135 <input type="checkbox"/>	K542 \$135 <input type="checkbox"/>	
<b>Additional battery</b> (recommended for institutional use)		K307 \$135 <input type="checkbox"/>	K307 \$135 <input type="checkbox"/>		
<b>Additional thigh straps</b>	CHOOSE	Narrow single	K318 \$50 <input type="checkbox"/>	K318 \$50 <input type="checkbox"/>	
		Narrow pair	K308 \$100 <input type="checkbox"/>	K308 \$100 <input type="checkbox"/>	
		Wide single	K328 \$50 <input type="checkbox"/>	K328 \$50 <input type="checkbox"/>	
		Wide pair	K338 \$100 <input type="checkbox"/>	K338 \$100 <input type="checkbox"/>	
<b>Swivel locks</b> (pair)		K319 \$185 <input type="checkbox"/>	K324 \$185 <input type="checkbox"/>		
<b>Removable back belt</b>		K327 \$60 <input type="checkbox"/>	K327 \$60 <input type="checkbox"/>		
<b>Total</b>		<b>K310</b> \$ _____	<b>K320</b> \$ _____		

\* Required


# TRAM dimensions

User dimensions (inches)		K310 and K320	
Girth		22 min – 60 max	
 <p><b>Key user dimension: girth</b></p> <p><b>Important:</b> User's weight or height must not exceed 350 lb or 76 inches</p>			
Item dimensions (inches)	K310 TRAM	K320 TRAM	
Overall length	45½	Same	
Overall width	27½–40	Same	
Overall height	43½–58	Same	
TRAM weight (lb)	70	Same	
Turning diameter	50	Same	
Floor to top of base leg for furniture clearance	6¾	4¼	
Min. user armpit height (when standing or ambulating)	30½	Same	
Max. user height (when standing or ambulating)	76	Same	
Max. working load (lb)	350	Same	

## Pelvic support

User dimensions (inches)			
Indicated by the red arrow below			
Key user dimension: half of thigh circumference measured at the groin.			
K306 small	K316 med.	K326 large	
6–10	8–14	10–18	
			
Product dimensions (inches)			
K306 small	K316 med.	K326 large	
<b>a</b>	8½	12	13½
<b>b</b>	16	20	25
<b>c</b>	9	12	14

## Hip positioner

Product dimensions (inches) – without pad		
		
	K667 small	K668 large
<b>a</b> Crotch length	2½	3½
<b>b</b> Width	1¼	1¾
<b>c</b> Length	9	11