

2022 order form


Use dimension chart to select appropriate size.




		K310 Standard (better for gait training)		K320 Low-base (fits under low beds)	
* Frame includes one battery, one battery charger, accessories tote and charger power cord. (International power cords available on request.)	CHOOSE	Frame with scale	K301 \$6,100 <input type="checkbox"/>	K321 \$6,100 <input type="checkbox"/>	
		Frame	K302 \$4,990 <input type="checkbox"/>	K322 \$4,990 <input type="checkbox"/>	
* Thigh straps (pair)	CHOOSE	Narrow	K308 \$160 <input type="checkbox"/>	K308 \$160 <input type="checkbox"/>	
		Wide	K338 \$160 <input type="checkbox"/>	K338 \$160 <input type="checkbox"/>	
Switch pole			K666 \$235 <input type="checkbox"/>	K666 \$235 <input type="checkbox"/>	
Arm supports (pair)	CHOOSE	Small arm prompts	K305 \$495 <input type="checkbox"/>	K305 \$495 <input type="checkbox"/>	
		Large arm prompts	K315 \$570 <input type="checkbox"/>	K315 \$570 <input type="checkbox"/>	
		Arm platforms	K397 \$370 <input type="checkbox"/>	K397 \$370 <input type="checkbox"/>	
		Arm platforms w/ handgrips	K396 \$420 <input type="checkbox"/>	K396 \$420 <input type="checkbox"/>	
Pelvic support	CHOOSE	Small	K306 \$250 <input type="checkbox"/>	K306 \$250 <input type="checkbox"/>	
		Medium	K316 \$285 <input type="checkbox"/>	K316 \$285 <input type="checkbox"/>	
		Large	K326 \$300 <input type="checkbox"/>	K326 \$300 <input type="checkbox"/>	
Hip positioner	CHOOSE	Small	K667 \$310 <input type="checkbox"/>	K667 \$310 <input type="checkbox"/>	
		Large	K668 \$340 <input type="checkbox"/>	K668 \$340 <input type="checkbox"/>	
Hip positioner pad	CHOOSE	Small	K541 \$105 <input type="checkbox"/>	K541 \$105 <input type="checkbox"/>	
		Large	K542 \$142 <input type="checkbox"/>	K542 \$142 <input type="checkbox"/>	
Additional battery (recommended for institutional use)			K307 \$145 <input type="checkbox"/>	K307 \$145 <input type="checkbox"/>	
Additional thigh straps	CHOOSE	Narrow single	K318 \$80 <input type="checkbox"/>	K318 \$80 <input type="checkbox"/>	
		Narrow pair	K308 \$160 <input type="checkbox"/>	K308 \$160 <input type="checkbox"/>	
		Wide single	K328 \$80 <input type="checkbox"/>	K328 \$80 <input type="checkbox"/>	
		Wide pair	K338 \$160 <input type="checkbox"/>	K338 \$160 <input type="checkbox"/>	
Swivel locks (pair)			K319 \$190 <input type="checkbox"/>	K324 \$190 <input type="checkbox"/>	
Total			K310 \$ _____	K320 \$ _____	

* Required


TRAM dimensions

User dimensions (inches)		K310 and K320	
Girth		22 min – 55 max	
 <p>Key user dimension: girth</p> <p>Important: User's weight or height must not exceed 350 lb or 76 inches</p>			
Item dimensions (inches)	K310 TRAM	K320 TRAM	
Overall length	45½	Same	
Overall width	27½–40	Same	
Overall height	43½–58	Same	
TRAM weight (lb)	70	Same	
Turning diameter	50	Same	
Floor to top of base leg for furniture clearance	6¾	4¼	
Min. user armpit height (when standing or ambulating)	30½	Same	
Max. user height (when standing or ambulating)	76	Same	
Max. working load (lb)	350	Same	

Pelvic support

User dimensions (inches) Indicated by the red arrow below			
Key user dimension: half of thigh circumference measured at the groin.			
K306 small	K316 med.	K326 large	
6–10	8–14	10–18	
			
Product dimensions (inches)			
K306 small	K316 med.	K326 large	
a	8½	12	13½
b	16	20	25
c	9	12	14

Hip positioner

Product dimensions (inches) – without pad		
		
	K667 small	K668 large
a Crotch length	2½	3½
b Width	1¼	1¾
c Length	9	11