

# 2019 order form

Use dimension chart to select appropriate size.

## Conversion kits

The tool-free conversion kit converts your small chair to medium as your child grows.

### Fabric without calf rest

Blue Z268 \$325  
 Pink Z269 \$325

(Includes seat, back, cover and chest strap)



### Fabric with calf rest

Blue Z266 \$400  
 Pink Z267 \$400


(Includes seat, back, cover and chest strap)



		Z210 small	Z220 medium	Z230 large	
<b>* Chair</b>		Z211 \$270 <input type="checkbox"/>	Z221 \$280 <input type="checkbox"/>	Z231 \$345 <input type="checkbox"/>	
<b>* Fabric cover</b>	CHOOSE	Blue	Z244 \$195 <input type="checkbox"/>	Z246 \$205 <input type="checkbox"/>	Z248 \$215 <input type="checkbox"/>
		Pink	Z245 \$195 <input type="checkbox"/>	Z247 \$205 <input type="checkbox"/>	Z249 \$215 <input type="checkbox"/>
<b>* Chest strap</b>	CHOOSE	With lateral positioning	Z212 \$75 <input type="checkbox"/>	Z222 \$85 <input type="checkbox"/>	Z232 \$95 <input type="checkbox"/>
		Without lateral positioning	Z213 \$65 <input type="checkbox"/>	Z223 \$65 <input type="checkbox"/>	Z233 \$85 <input type="checkbox"/>
<b>Head blocks</b>		Z214 \$92 <input type="checkbox"/>	Z214 \$92 <input type="checkbox"/>	Z234 \$92 <input type="checkbox"/>	
<b>Leg straps</b>		Z215 \$65 <input type="checkbox"/>	Z225 \$75 <input type="checkbox"/>	Z235 \$85 <input type="checkbox"/>	
<b>Calf rest</b> (calf rest must be chosen at the time of the original purchase. It cannot be added later.)		Z216 \$120 <input type="checkbox"/>	Z226 \$135 <input type="checkbox"/>	Z236 \$150 <input type="checkbox"/>	
<b>Tub stand</b>		Z217 \$270 <input type="checkbox"/>	Z217 \$270 <input type="checkbox"/>	Z217 \$270 <input type="checkbox"/>	
<b>Shower stand</b>		Z218 \$730 <input type="checkbox"/>	Z218 \$730 <input type="checkbox"/>	Z218 \$730 <input type="checkbox"/>	
<b>Tub transfer base</b> (calf rest recommended)	CHOOSE	Counterclockwise into tub	Z227 \$860 <input type="checkbox"/>	Z227 \$860 <input type="checkbox"/>	Z227 \$860 <input type="checkbox"/>
		Clockwise into tub	Z228 \$860 <input type="checkbox"/>	Z228 \$860 <input type="checkbox"/>	Z228 \$860 <input type="checkbox"/>
<b>Legs for tub transfer base</b>	CHOOSE	Standard	Z224 \$85 <input type="checkbox"/>	Z224 \$85 <input type="checkbox"/>	Z224 \$85 <input type="checkbox"/>
		Long (tub edge over 20" high)	Z238 \$85 <input type="checkbox"/>	Z238 \$85 <input type="checkbox"/>	Z238 \$85 <input type="checkbox"/>
<b>Additional head blocks</b>		Z214 \$92 <input type="checkbox"/>	Z214 \$92 <input type="checkbox"/>	Z234 \$92 <input type="checkbox"/>	
<b>Additional fabric cover</b>		Z219 \$195 <input type="checkbox"/>	Z229 \$205 <input type="checkbox"/>	Z239 \$215 <input type="checkbox"/>	
<b>Additional chest strap</b>	CHOOSE	With lateral positioning	Z212 \$75 <input type="checkbox"/>	Z222 \$85 <input type="checkbox"/>	Z232 \$95 <input type="checkbox"/>
		Without lateral positioning	Z213 \$65 <input type="checkbox"/>	Z223 \$65 <input type="checkbox"/>	Z233 \$85 <input type="checkbox"/>
<b>Total</b>		<b>Z210</b> \$ _____	<b>Z220</b> \$ _____	<b>Z230</b> \$ _____	

\* Required

# Wave dimensions

User dimensions (inches)	Z210 small	Z220 medium	Z230 large
Height	up to 46	40–56	50–74
 <p><b>Key user dimension: height</b> Select the appropriate bath chair by the user's overall height. Choose the size that allows for growth.</p> <p><b>Important:</b> Make sure that seat depth, width and back height are adequate. The user's weight must not exceed the maximum working load.</p>			
	Item dimensions (inches)	Z210 small	Z220 medium
Seat depth	10½	13	17
Seat width	15	15	17½
Back height	23	26½	32
Seat height	2–7	2–7	2–7
Seat height (incl. tub stand)	12½–18½	12½–18½	12½–18½
Seat height (incl. shower stand)	27½	27½	27½
Seat height (incl. tub transfer base) with long leg inserts	21–26 24½–29½	21–26 24½–29½	21–26 24½–29½
Length of mesh area when flat (without calf rest)	35	41	50½
Length of mesh area when flat (with calf rest)	42½	49	61
Base width	18	18	21
Base length (legs down) with legs flat	24 29½	24 29½	24 29½
Tub stand height (legs flat)	8	8	8
Tub stand length (legs flat)	39	39	39
Tub stand width (legs flat)	17½	17½	17½
Shower stand height	22	22	22
Shower stand width	27	27	27
Shower stand length	34	34	34
Tub transfer base footprint	28 x 30	28 x 30	28 x 30
Max tub height for transfer base with long leg inserts	20 24	20 24	20 24
Maximum working load (lbs)	75	150	250
Back and seat angle adjustments		0°, 15°, 30°, 45°, 60°, 75°	



All three sizes are available in blue and pink.