


# 2019 order form

Use dimension chart to select appropriate size.



User dimensions (inches)	E420 small	E430 large
Height	30–50	46–72
 <b>Key user dimension: height</b> Select the appropriate stander by the user's overall height. Choose the model that allows for growth.		
	<b>Important:</b> User's weight must not exceed the maximum working load.	
Item dimensions (inches)	E420 small	E430 large
Board length and width	53 x 16½	73 x 19
Distance between laterals	7–13	9–16
Height when horizontal	28½	21–30
Height when vertical	54½	74
Base length and width	37½ x 27	46½ x 29½
Max. working load (lbs)	100	250

\* **Supine Stander** includes padding, footboard, adjustable armrests, hip strap, chest strap, leg straps, one pair of laterals, one pair of head blocks and double-locking casters.

		E420 small		E430 large	
		E40	\$2,705 <input type="checkbox"/>	E41	\$4,265 <input type="checkbox"/>
<b>Additional laterals (pair)</b>	<b>CHOOSE</b> Small	E504	\$140 <input type="checkbox"/>	N/A	
	Large	N/A		E439	\$205 <input type="checkbox"/>
<b>Tray</b>		E404	\$270 <input type="checkbox"/>	E414	\$335 <input type="checkbox"/>
<b>CHOOSE</b>	<b>Adjustable abduction wedge</b>	E945	\$225 <input type="checkbox"/>	E945	\$225 <input type="checkbox"/>
	<b>Round abductor</b>	E933	\$95 <input type="checkbox"/>	E933	\$95 <input type="checkbox"/>
<b>Collar</b>		E934	\$75 <input type="checkbox"/>	E934	\$75 <input type="checkbox"/>
<b>Sandals (pair)</b>	<b>CHOOSE</b> Small	E422	\$85 <input type="checkbox"/>	N/A	
	Medium	E433	\$94 <input type="checkbox"/>	E433	\$94 <input type="checkbox"/>
	Large	N/A		E432	\$97 <input type="checkbox"/>
<b>Wedges (pair, requires sandals)</b>	<b>CHOOSE</b> Small	E424	\$61 <input type="checkbox"/>	N/A	
	Medium	E435	\$67 <input type="checkbox"/>	E435	\$67 <input type="checkbox"/>
	Large	N/A		E434	\$70 <input type="checkbox"/>
<b>Hand Anchor</b> (For more hand/arm support options see pp. 98–101.)		K821	\$145 <input type="checkbox"/>	K821	\$145 <input type="checkbox"/>
<b>Total</b>		<b>E420</b>	<b>\$ _____</b>	<b>E430</b>	<b>\$ _____</b>

\* Required