


2021 order form

Use dimension chart to select appropriate size.



User dimensions (inches)	E420 small	E430 large
Height	30–50	46–72
 Key user dimension: height Select the appropriate stander by the user's overall height. Choose the model that allows for growth.		
	Important: User's weight must not exceed the maximum working load.	
Item dimensions (inches)	E420 small	E430 large
Board length and width	53 x 16½	73 x 19
Distance between laterals	7–13	9–16
Height when horizontal	28½	21–30
Height when vertical	54½	74
Base length and width	37½ x 27	46½ x 29½
Max. working load (lb)	100	250

* **Supine Stander** includes padding, footboard, adjustable armrests, hip strap, chest strap, leg straps, one pair of laterals, one pair of head blocks and double-locking casters.

		E420 small	E430 large
		E40	\$2,840 <input type="checkbox"/>
		E41	\$4,450 <input type="checkbox"/>
Additional laterals (pair)	CHOOSE		
	Small	E504	\$145 <input type="checkbox"/>
	Large	N/A	E439 \$215 <input type="checkbox"/>
Tray		E404	\$285 <input type="checkbox"/>
	E414	\$350 <input type="checkbox"/>	
CHOOSE	Adjustable abduction wedge	E945	\$235 <input type="checkbox"/>
	Round abductor	E933	\$100 <input type="checkbox"/>
	E933	\$100 <input type="checkbox"/>	
Collar		E934	\$80 <input type="checkbox"/>
	E934	\$80 <input type="checkbox"/>	
Sandals (pair)	CHOOSE		
	Small	E422	\$90 <input type="checkbox"/>
	Medium	E433	\$100 <input type="checkbox"/>
	Large	N/A	E432 \$100 <input type="checkbox"/>
Wedges (pair, requires sandals)	CHOOSE		
	Small	E424	\$65 <input type="checkbox"/>
	Medium	E435	\$70 <input type="checkbox"/>
	Large	N/A	E434 \$75 <input type="checkbox"/>
Hand Anchor (For more hand/arm support options see pp. 104–107.)		K821	\$150 <input type="checkbox"/>
	K821	\$150 <input type="checkbox"/>	
Total		E420	\$ _____
		E430	\$ _____

* Required