


# Support Station dimensions

User dimensions (inches)		K710
Height		46 – 80
	<b>Important:</b> User height and weight must not exceed the maximum.	
Item dimensions (inches)		K710
Trunk board size		24 x 24
Maximum trunk board height		45
Minimum trunk board height		26
Height of frame		65½
Width of frame		27½
Depth when trunk board horizontal		31
Depth when collapsed against wall		8
Max. working load (lbs)		250

## 2019 order form

Use dimension chart to select appropriate size.



\* Frame

		K710	
* Frame		K715	\$3,120 <input type="checkbox"/>
Support strap		K716	\$95 <input type="checkbox"/>
Kneeboard		K717	\$615 <input type="checkbox"/>
Handholds	Pair	K718	\$76 <input type="checkbox"/>
	Single	K719	\$38 <input type="checkbox"/>
<b>Total</b>		<b>K710</b>	<b>\$ _____</b>

\* Required