


Support Station dimensions

User dimensions (inches)		K710
Height		46 – 80
 <p>Important: User height and weight must not exceed the maximum number.</p>		
Item dimensions (inches)		K710
Trunk board size		24 x 24
Maximum trunk board height		45
Minimum trunk board height		26
Height of frame		65½
Width of frame		27½
Depth when trunk board horizontal		31
Depth when collapsed against wall		8
Max. working load (lb)		250

2021 order form



* Frame

		K710	
		K715	\$3,270 <input type="checkbox"/>
Support strap		K716	\$100 <input type="checkbox"/>
Kneeboard		K717	\$645 <input type="checkbox"/>
Handholds	Pair	K718	\$80 <input type="checkbox"/>
	Single	K719	\$40 <input type="checkbox"/>
Total		K710	\$ _____

* Required