


# 2019 order form

Use dimension chart to select appropriate size.



User dimensions (inches)	R120 small	R130 medium	R140 large
Inside leg length	17 – 22	22 – 28	25 – 35
 <p><b>Key user dimension: inside leg length</b> The extended leg should reach from seat to pedal comfortably when both feet and torso are strapped into place.</p> <p><b>Important:</b> User's weight must not exceed the maximum working load. Rifton tricycles are not intended for clients with strong reflexes or poorly developed trunk balance and head control. Consult your client's therapist before ordering.</p>			
	Item dimensions (inches)	R120 small	R130 medium
Overall height (at lowest position)	28½	31	36½
Overall length	46	52	65
Overall width	27½	27½	31½
Wheel diameter	13	15	19
Seat center to extended pedal	17 – 22	22 – 28	25 – 35
Seat to top of trunk support	15½ – 22	15½ – 22	15½ – 22
Gear ratio (approximate)	1:1	1:1	1:1
Item weight (lbs)	47	56	65
Distance between trunk laterals	8 – 14	8 – 14	8 – 14
Max. working load (lbs)	150	160	200

		R120 small	R130 medium	R140 large	
<p>* <b>Tricycle</b> frame includes self-leveling pedals with straps, backrest tube with handle, tote box for storage and handbrake for the large only.</p> <p><b>Specify color:</b> Red, blue, pink or lime.</p>	CHOOSE	R112     \$1,010 <input type="checkbox"/>	R113     \$1,150 <input type="checkbox"/>	R114     \$1,945 <input type="checkbox"/>	
		Red <input type="checkbox"/>	Red <input type="checkbox"/>	Red <input type="checkbox"/>	
		Blue <input type="checkbox"/>	Blue <input type="checkbox"/>	Blue <input type="checkbox"/>	
		Pink <input type="checkbox"/>	Pink <input type="checkbox"/>	Pink <input type="checkbox"/>	
		Lime <input type="checkbox"/>	Lime <input type="checkbox"/>	Lime <input type="checkbox"/>	
<p>* <b>Seat</b> (includes seatbelt)</p>	CHOOSE	Small	R122     \$205 <input type="checkbox"/>	R132     \$205 <input type="checkbox"/>	N/A
	Large	R128     \$350 <input type="checkbox"/>	R138     \$350 <input type="checkbox"/>	(included)	
<p>* <b>Handlebar</b></p>	CHOOSE	Straight	R126     \$130 <input type="checkbox"/>	R136     \$140 <input type="checkbox"/>	R146     \$220 <input type="checkbox"/>
	Loop	R129     \$130 <input type="checkbox"/>	R139     \$140 <input type="checkbox"/>	R149     \$220 <input type="checkbox"/>	
<p>* <b>Back support</b></p>	CHOOSE	Trunk support system	R147     \$475 <input type="checkbox"/>	R147     \$475 <input type="checkbox"/>	R147     \$475 <input type="checkbox"/>
	Backrest pad	R145     \$105 <input type="checkbox"/>	R145     \$105 <input type="checkbox"/>	R145     \$105 <input type="checkbox"/>	
<p><b>Headrest</b> (requires trunk support system)</p>	CHOOSE	Winged	R159     \$225 <input type="checkbox"/>	R159     \$225 <input type="checkbox"/>	R159     \$225 <input type="checkbox"/>
	Contoured	R156     \$215 <input type="checkbox"/>	R156     \$215 <input type="checkbox"/>	R156     \$215 <input type="checkbox"/>	
	Flat	R157     \$205 <input type="checkbox"/>	R157     \$205 <input type="checkbox"/>	R157     \$205 <input type="checkbox"/>	
<p><b>Handbrakes</b> (Must be factory installed. Cannot be added later.)</p>		N/A	R125     \$105 <input type="checkbox"/>	(included)	
<p><b>Communication tray</b></p>	CHOOSE	Small	K225     \$185 <input type="checkbox"/>	K225     \$185 <input type="checkbox"/>	K225     \$185 <input type="checkbox"/>
	Large	N/A	K245     \$190 <input type="checkbox"/>	K245     \$190 <input type="checkbox"/>	
<b>Abductor</b>		R121     \$90 <input type="checkbox"/>	R131     \$145 <input type="checkbox"/>	R141     \$165 <input type="checkbox"/>	
<b>Front guide bar</b>		R135     \$130 <input type="checkbox"/>	R135     \$130 <input type="checkbox"/>	R135     \$130 <input type="checkbox"/>	
<b>Rear steering bar</b>		R123     \$230 <input type="checkbox"/>	R123     \$230 <input type="checkbox"/>	R143     \$255 <input type="checkbox"/>	
<b>Front pulley</b>		R144     \$43 <input type="checkbox"/>	R144     \$43 <input type="checkbox"/>	R144     \$43 <input type="checkbox"/>	
<b>Additional backrest pad</b>		R145     \$105 <input type="checkbox"/>	R145     \$105 <input type="checkbox"/>	R145     \$105 <input type="checkbox"/>	
<b>Stationary stand</b>		R155     \$115 <input type="checkbox"/>	R155     \$115 <input type="checkbox"/>	R155     \$115 <input type="checkbox"/>	
<b>Total</b>		<b>R120</b> \$ _____	<b>R130</b> \$ _____	<b>R140</b> \$ _____	

\* Required